 ceramill map400 **CTOMAPT**
DENTAL SOLUTIONS**Betriebsanleitung****User Manual****Mode d'emploi****Istruzioni d'uso****Modo de empleo**

Deutsch	3 - 27
English	28 - 51
Français	52 - 77
Italiano	78 - 102
Español	103 - 128

**AMANN GIRRBACH**

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Tab. 5

- Translation of the original Instruction Manual -

Table of Contents

1	Explanation of Symbols	29	6	Application Tips	49
			6.1	Fixator	49
			6.2	Positioning the Model	49
2	General Safety Instructions	29			
3	Suitable Personnel	29	7	Registration and Download Information	50
			7.1	Registration, M-center	50
			7.2	Download Infos	50
4	Machine Specifications	30	8	Cleaning and Maintenance	50
4.1	Delivery Scope of the Complete System (incl. Virtual Articulator) .	30	8.1	Cleaning	50
4.2	Intended Use	30	8.2	Maintenance	50
4.3	CE Declaration of Conformity	30			
4.4	Components and Interfaces	31	9	Malfunctions, Repairs and Warranty ..	50
4.4.1	Ceramill Map400	31	9.1	Malfunctions	50
4.4.2	Ceramill Transferkit	31	9.2	Repairs	50
4.4.3	PC Model T3500	32	9.3	Warranty	50
5	Installation	33	10	Environmental Protection	50
5.1	Transport	33			
5.2	Setting Up	33	11	Technical Data	51
5.3	Installation of the PC	33			
5.4	Installation of the Software on the PC	33			
5.4.1	Installation of the Anti-virus Program	33			
5.4.2	Installation of the Ceramill Mind and Ceramill Map Software	34			
5.5	Installation of the Scanner	34			
5.6	Activating Additional Modules	34			
5.7	Calibration of Axes	35			
5.8	Calibrating the Artex Fixator	35			
5.9	Registration of the Artex Fixator .	36			

1 Explanation of Symbols

Warning indications

Warning indications in the text are marked with a colour-backed triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- **NOTE** means that property damage can occur.
- **CAUTION** means that light to fairly serious personal injury can occur.
- **WARNING** means that serious personal injury can occur.
- **DANGER** means that serious personal injury can occur.

Important information

Important information that do not lead to hazards for humans or property damage are marked with the icon aside and are also boxed.

Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
–	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

Tab. 1

Other symbols on the machine

Symbol	Meaning
	USB port (socket type B)
	Fuse

Tab. 2

2 General Safety Instructions

When installing, starting-up and operating the machine, always observe the following safety instructions:

CAUTION:

A defective machine can lead to malfunctions!
On detection of damage or a functional defect of the machine:

- ▷ Label the machine as defective.
- ▷ Prevent further operation until the machine has been repaired.

NOTE:

▷ Switch the machine off when not in use or unsupervised for longer periods, e.g., overnight. This method of saving electrical energy also benefits the environment.

3 Suitable Personnel

NOTE:

Starting-up and operation of the machine may only be carried out by trained specialised personnel.

4 Machine Specifications

4.1 Delivery Scope of the Complete System (incl. Virtual Articulator)

- _ Scanner Ceramill Map400 (179140)
 - Transferkit
 - Measuring template
 - USB cable (2 x)
 - Mains cable (1 x)
 - Operating Instructions
 - CD with installation instructions and calibration data
 - Calibration model
 - _ PC (179170)
 - Windows 7 operating system
 - Mains cable (1 x)
 - Monitor (1 x VGA cable, 1 x DVI cable, 1 x mains cable)
 - Keyboard, mouse, mousepad, headset
 - _ Ceramill Mind Package (179150)
 - Ceramill Mind Dongle
 - Quickguide Installation
 - Ceramill Mind CD
 - DVD with application tutorials
- ▷ After unpacking, check the machine for completeness and possible transport damages. Please claim any transport damages immediately with your supplier.

4.2 Intended Use

The Ceramill Map400 is a PC-controlled 3D scanner for scanning dental tooth models.

NOTE:

The scan result depends on the materials being scanned. Metal surfaces, for example, can lead to artefacts. Scan spray can improve the scan capability of materials.

For improper use of the system, AmannGirrbach shall not assume any liability whatsoever.

Unauthorised modifications/alterations of the machine shall also void the warranty.

4.3 CE Declaration of Conformity

In terms of design and performance, this product complies with the European Directives and the supplementary national requirements. Conformity has been confirmed with the CE marking.

The product's Declaration of Conformity can be viewed under www.amanngirrbach.com.

4.4 Components and Interfaces

4.4.1 Ceramill Map400


Fig. 1 Machine overview, front

- 1 Artex Scan base plate
- 2 Operation button


Fig. 2 Machine overview, rear

- 1 USB port (control system)
- 2 Main switch
- 3 Functional earthing
- 4 USB port (camera)
- 5 Power supply connection

4.4.2 Ceramill Transferkit


Fig. 3

- A Scan set-up for models without articulation base
- B Scan set-up for models with articulation base
- C Scan set-up for two models in relation to the articulator
- 1 No plast model retainer
- 2 Model plate
- 3 Screws with various lengths, depending on application of spacer plates
- 4 Allen key
- 5 Fixation plate with clamp
- 6 Spacer plate 2 (20 mm); optional, depending on model height
- 7 Spacer plate 1 (10 mm); optional, depending on model height
- 8 Artex fixator

The Artex scan base plate is stationary in the Ceramill Map400 scanner. The spacer plates **1** **7** and **2** **6** and the fixation plate with clamp **5** are screwed to the Artex scan base plate as required. The choice of spacer plates and screw **3** depend on the height of the model to be scanned.

4.4.3 PC Model T3500


Fig. 4

- 1** Power supply connection, PC
- 2** Network connection, Motion
- 3** Digital connection, monitor
- 4** Network connection, Internet
- 5** USB connections for keyboard, mouse, scanner and dongle

5 Installation

5.1 Transport

- ▷ For transport, have two persons lift the equipment by the floor plate.

5.2 Setting Up

- The machine is intended exclusively for use within dry, closed rooms.
- Minimum space requirements for the Ceramill Map400 (incl. connections):
 - Width: 420 mm
 - Depth: 450 mm (excl. connection cables)
 - Height: 610 mm
- The set-up surface must withstand a load of 35 kg.
- During operation, the room temperature must be between 15 °C and 30 °C; avoid large variations in temperature.
- ▷ Set up the machine on a sturdy work table or work bench (off the ground and clear of walls; important due to natural oscillation).
- ▷ The opening of the scan chamber must face away from windows and strong artificial lighting.
- ▷ Adjust the height of the machine by screwing the feet in or out, ensuring that the machine is positioned firmly on all four feet.

5.3 Installation of the PC

- ▷ Connect the mains cable to the power supply connection of the PC as well as to a mains socket outlet.
- ▷ Connect a network cable to the network connection of the PC and the Internet router/modem.
- ▷ Unplug the blue VGA cable (analog) from the monitor.

- ▷ Connect the white DVI cable (digital) to the monitor as well as to the digital connection of the PC.
- ▷ Connect the mouse and the keyboard to the PC.

5.4 Installation of the Software on the PC

- ▷ Switch the scanner on via the main switch.
- ▷ Switch the PC on.

The installation shall take place under the special user "Administrator":

– **Username: Administrator**

– **Password: ceramill mall**

- ▷ Use this user only for the installation.

- ▷ For the installation, change the user and login as Administrator.

After restarting, the PC logs "cadcam" as the default user. A password is not factory-assigned.

Language Selection

The factory default language for the Ceramill Digital System is English. For selection of the desired language after booting the PC:

- ▷ On the desktop, open *Manual Language Selector* and follow the instructions.

5.4.1 Installation of the Anti-virus Program

For protection of the PC:

- ▷ On the desktop, click on *Microsoft Security Essentials*.
This link takes you to the Microsoft Security Essentials homepage.
- ▷ Follow the instructions on the homepage to install the anti-virus program.

5.4.2 Installation of the Ceramill Mind and Ceramill Map Software

The PCs are factory-shipped with a basic installation. The Ceramill Mind and Ceramill Map software are to be installed by the user in this sequence. For the installation, login as Administrator is required:

- ▷ Insert the Ceramill Mind CD provided.
- ▷ Open the installation instructions on the CD and follow the instructions.
- ▷ Insert the Ceramill Map CD provided.
- ▷ Open the installation instructions on the CD and follow the instructions.

To switch to the desired language for the installed software components:

- ▷ On the desktop, open *Manual Language Selector*, follow the instructions and select the installed scanner.

5.5 Installation of the Scanner

As each scanner is individually factory-calibrated, the data on the installation CD are only valid for the respective scanner.

Prior to installation:

- ▷ Check if the serial number of the scanner on the type plate is identical with the serial number stated on the CD!

Prior to installation:

- ▷ Make sure that the main switch of the scanner is positioned on **O**.
- ▷ Connect the two enclosed USB cables to the scanner and the PC.
- ▷ Plug the Ceramill Mind Dongle in a USB port on the rear side of the PC.
- ▷ Connect the mains cable to the power supply connection of the scanner as well as to a mains socket outlet.

5.6 Activating Additional Modules

Additional modules must be activated in the scanner software before they can be used.

- ▷ Click on *Options > Access > Activation Code ...*, enter the received scanner-specific activation code.

- ▷ Restart the software to activate the module.

5.7 Calibration of Axes

Following the set-up and installation of the Ceramill Map400 scanner to the PC, the scanner must be calibrated using the supplied calibration model.

- Mount spacer plate No. 1, the fixation plate with the clamp and the model plate (please refer to section 4.4.2).

Fig. 5 Calibration model

- Insert the calibration model
- Create a new patient case and start Ceramill Map from the data base.
- Cancel the subsequent command prompt for the 2D scan.
- Start the axes calibration under *Options > Service > Axis calibration*.
- Follow the software instructions.

If required:

- Re-enter the calibration-model values entered during installation under *Options > Service > Register calibration object....*

The scanner should be recalibrated at regular intervals.

5.8 Calibrating the Artex Fixator

The calibration of the fixator with the Splitex calibration key should be checked at regular intervals.

Before the Artex Fixator is used in the scanner, it must be calibrated with the Splitex calibration key:


Fig. 6

- 1 Fixator, bottom element
- 2 Adhesive surface
- 3 Splitex plate
- 4 Splitex calibration key
- 5 Fixator, upper element
- 6 Lateral screw
- 7 Rear screw

- Open the Fixator.
- Place Splitex calibration key **4** onto upper element **5**.
- Place Splitex plate **3** onto calibration key.
- Apply adhesive (Loctite) to bottom side of the Splitex plate.
- Spray activator onto bottom element of fixator **1**.

- ▷ Tilt fixator together to the stop and lock it with the lateral screw **6**.
- ▷ Allow adhesive connection to cure according to adhesive instructions.
- ▷ Open fixator.
- ▷ Remove the Splitex calibration key.

5.9 Registration of the Artex Fixator

When using the fixator for the first time, the upper and lower jaw scans must be assigned manually to the vestibular scan by means of the 3-point method.

This manual assigning is used for the latter automatic pre-alignment of the jaw scans. The manual method can be repeated any time, as required.

The scanner must be calibrated before using for the first time (please refer to section 5.7)!

NOTE:

Faulty results!

- ▷ Check the articulator being used to articulate the models with the Splitex calibration key and recalibrate, as required (see "Calibration Splitex" Manual).
- ▷ For articulation of the jaw models used for registration, set the supporting pin of the articulator to zero.

Requirements of the models used for the initial 3-point assignment:

- _ The models must have been plastered into the Artex CR articulator.
- _ Sufficient interlocking and supporting of upper and lower jaw to each other, to ensure successful allocation.

NOTE:

Faulty registration of the fixator!

- ▷ Always place jaw models with plaster base in the scanner. Do not use models with Artex Noplast.

- ▷ Set three model marks using a small router (recommended).


Fig. 7

- ▷ Start the Ceramill Mind data base.
- ▷ Create a fictive patient case with indications in the upper and lower jaw (for example, coping in upper jaw and antagonist in lower jaw).
- ▷ Select *Artex CR* as the antagonist type.


Fig. 8

- ▷ Start the scanner.
- ▷ Using the required spacer plate(s) (0, 1, 2 or 1+2), position the upper jaw within the measuring range in the scanner (please refer to section 6.2).
- ▷ Perform a 2D scan.

- ▷ Select the number of spacer plates.
- ▷ Under Scan type, select *Complete Jaw*.


Fig. 9

- ▷ Start the 3D scan with the *Next* button.

▷ Check the scan result and re-scan, if required.


Fig. 10

- ▷ Follow the software instructions.
- ▷ Remove the upper jaw.
- ▷ Using the required spacer plate (0, 1, 2 or 1+2), position the lower jaw within the measuring range in the scanner (please refer to section 6.2).
- ▷ Perform a 2D scan.
- ▷ Select the number of spacer plates.
- ▷ Under Scan type, select *Complete Jaw*.
- ▷ Start the 3D scan with the *Next* button.
- ▷ Check the scan result and re-scan, if required.
- ▷ Follow the software instructions.
- ▷ Remove the lower jaw together with the model set-up (see figure 3 on page 31).
- ▷ Position both models in the fixator.

NOTE:

The tiltable upper element of the fixator may not be put into the scanner when not secured (affixed)!

▷ Firmly tighten the screw on the rear side of the fixator!

▷ Afterwards, tighten the lateral screw.

▷ Place the fixator into the scanner.

The clamp of the fixation plate must engage in the bottom element of the fixator.


Fig. 11 Fixator in the Ceramill Map400

1 Artex fixator

2 Screw


3 Fixation plate with clamp

4 Spacer plate (optional)

5 Artex Scan base plate

- Start the vestibular scan by clicking on the *Next* button.
- Follow the software instructions.
- The software requires a 3-point assignment for the necessary registration of the fixator.
- Confirm with *OK*.
- Define at least three corresponding points for registration on the lower-jaw and vestibular scan.

Objective	Action
Setting a point	Press the left mouse button 
Turning the model	Press the <i>Strg/Ctrl</i> key and the right mouse button and move the mouse 
Moving the model	Press the <i>Strg/Ctrl</i> key and both mouse buttons and move the mouse 

Tab. 3

► Set three points in the lower jaw (on the model marks).

The last set point can be deleted by pressing the *Entf/Del* key.


Fig. 12

> Switch to the vestibular scan by pressing the *space bar*.


Fig. 13

► There, set three points at the same marks in the same sequence.

The lower jaw is assigned.

In case the assignment was created incorrectly, it can be repeated with the *Articulation* function.


Fig. 14

> Start the upper jaw scan by clicking on the *Next* button.


Fig. 15

- Carry out the registration of the upper jaw in the same manner as with the lower jaw.
The upper jaw is assigned.


Fig. 16

► Start the fine alignment by clicking on the *Nexf* button.


Fig. 17

After fine alignment of both jaws:

▷ Start the matching process by clicking on the *Next* button.

The defined points for the fixator registration are stored in the software.


Fig. 18

The registration of the fixator is completed. From now on, all subsequent jaw scans of patient cases are automatically assigned.

The more accurate the points are set, the more precise the assignment. Small deviations are permitted, and will be corrected in the matching process!

When the assignment was created incorrectly:

- ▷ Select *Options > Service > Reset jaw alignment* to reset the faulty registration.
- ▷ Click on the *Articulation* button.
- ▷ Select *Manual alignment*
- ▷ Start the registration process by clicking on *OK*.

6 Application Tips

6.1 Fixator

When using the fixator for patient cases with sufficient remaining dentition (assured occlusion), thus ensuring sufficient jaw support:

- ▷ Do not tighten the rear screw of the fixator.

This ensures maximum, optimal interlocking of the jaws during the scan process.

The tiltable upper element of the fixator may not be put into the scanner when not secured (affixed)!

- ▷ Always tighten the lateral screw of the fixator.

When using the fixator for patient cases where no support of the jaws can be achieved through the remaining dentition (terminal occlusion):

- ▷ Articulate the models with the supporting pin position set to zero.
- ▷ Tighten the rear screw of the fixator.

NOTE:

Faulty automatic assignment of jaw models!

- ▷ When using the fixator (also applies for jaw scans), only place models with plaster base in the scanner!

6.2 Positioning the Model

The model must be positioned with the spacer plates in such a manner, that it is within the scan focus. In the scanner, the scan focus is located on the level of the tilting axis. For correct positioning, it is recommended to use the provided measuring template.

- ▷ Position the measuring template on the upper side of the tilting axis.
- ▷ Check if the areas of the jaw model to be scanned are within the marked measuring range.


Fig. 19 Using the measuring template

Using the No Plast Model Retainer

- ▷ Affix the model in the No Plast model retainer.
- ▷ Loosen the screws in the center part of the No Plast model retainer and align the model in the scan focus.
- ▷ Lock the position with the knurled screw and then additionally tighten the Allen screw with the Allen key.

7 Registration and Download Information

7.1 Registration, M-center

To download updates and upload data to the manufacturing center, it is necessary to register on the M-Center page. To register, log on to www.ceramill-m-center.com. Information and instructions for registration can also be found at the specified address.

7.2 Download Infos

Further instructions (video tutorials) on operating the Ceramill Map scanner and the Ceramill Mind software can be downloaded under www.ceramill-m-center.com.

Software updates are also available on our website under www.ceramill-m-center.com.

8 Cleaning and Maintenance

8.1 Cleaning

The scanner must be cleaned regularly.

- ▷ Switch the scanner off.

NOTE:

Possible damage to the 3D sensor:

- ▷ Do **not** clean the optics of the 3D sensor in the top of the scanner!

- ▷ Carefully vacuum out the interior of the scanner.
- ▷ Wipe off the housing with a moist cloth. Do not use aggressive cleaning agents.

8.2 Maintenance

Maintenance may be carried out only by trained, specialized personnel.

9 Malfunctions, Repairs and Warranty

9.1 Malfunctions

In case of malfunctions:

- ▷ Restart the software.
- ▷ Restart the scanner.
- ▷ Restart the PC.
- ▷ Check the cable connections; disconnect and then re-connect, as required.

If these measures do not correct the problem:

- ▷ Contact the Ceramill helpdesk (workdays from 8 to 5 p.m.).
 - Germany: +49 7231 957 100
 - Austria: +43 5523 62333 390
 - International: +43 5523 62333 399
 - helpdesk@amanngirrbach.com

9.2 Repairs

Repairs may be carried out only by trained, specialized personnel.

9.3 Warranty

The warranty complies with the statutory provisions. For more information, please refer to our General Terms of Business.

10 Environmental Protection

Packaging

In terms of packaging, AmannGirrbach participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally-friendly and recyclable.

Waste Equipment

Waste equipment contain materials that should be sorted for environmental-friendly recycling.

- After the service life has elapsed, dispose of the scanner via the public disposal systems in an environmentally-friendly manner.

11 Technical Data

Subject to changes.

Technical Data - Ceramill Map400

	Unit	Ceramill Map400
Art. No.	–	179140
Dimensions (D × W × H)		
_ closed	mm	415 × 407 × 535
_ open	mm	415 × 407 × 606
Weight	kg	31
Electrical connections		
_ Voltage	V/Hz	100-240/50-60
_ Fuse	A	2 × 1,25
Power output	W	80
Accuracy	µm	< 20
Allowable temperature range (operation)	°C	15 - 30
Interfaces	–	USB

Tab. 4

Technical Data of the PC

	Unit	Model T3500
Art. No.	–	179170
Weight	kg	26
Processor type/ cycle		Intel Xenon Quad/ 2.67GHz
Random access memory	GB RAM	12
Installed fixed-disk storage:	GB	232
Graphic card		Nvidia Quadro 2000
Installed operating system	–	Windows 7 Ultimate 64 bit, Recovery Software
Virus protection	–	Microsoft Security Essentials
Monitor:		
_ Size	Inch	22
_ Resolution	Pixel	1680 × 1050

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Amann Girrbach AG
Amann Girrbach GmbH