# Hygoclave 50 Hygoclave 50 Plus

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# Hygoclave 50 Hygoclave 50 Plus

EN

Installation and operating instructions

6045050000

# Contents

## Important information

|                                                            |           |
|------------------------------------------------------------|-----------|
| <b>1 Safety standards .....</b>                            | <b>43</b> |
| 1.1 Intended users.....                                    | 43        |
| 1.2 Intended use.....                                      | 43        |
| 1.3 General safety warnings.....                           | 43        |
| 1.4 Operating and transport environmental conditions ..... | 44        |
| 1.5 Environmental disposal instructions.....               | 44        |

## Product description

|                                                                    |           |
|--------------------------------------------------------------------|-----------|
| <b>2 Equipment description .....</b>                               | <b>45</b> |
| 2.1 Variants.....                                                  | 46        |
| 2.2 Main features .....                                            | 46        |
| 2.2.1 Sterilisation programs .....                                 | 46        |
| 2.2.2 Test programs .....                                          | 46        |
| 2.2.3 Water filling systems .....                                  | 46        |
| 2.2.4 Automatic sterilisation parameter control .....              | 46        |
| 2.2.5 Motorised closure .....                                      | 46        |
| 2.2.6 Bio-X filter.....                                            | 46        |
| 2.2.7 USB port .....                                               | 46        |
| 2.2.8 Serial port - RS232.....                                     | 46        |
| 2.2.9 Touch screen display .....                                   | 46        |
| 2.2.10 LAN connection.....                                         | 46        |
| 2.2.11 Safety devices .....                                        | 46        |
| 2.3 Accessories included in the medical device certification ..... | 47        |

## Installation

|                            |           |
|----------------------------|-----------|
| <b>3 Installation.....</b> | <b>48</b> |
| 3.1 First power-up.....    | 49        |

## Operating instructions

|                                                                |           |
|----------------------------------------------------------------|-----------|
| <b>4 Operating instructions.....</b>                           | <b>50</b> |
| 4.1 Symbols .....                                              | 51        |
| 4.2 Menu architecture.....                                     | 54        |
| 4.3 Selecting sterilisation programs.....                      | 55        |
| 4.3.1 Sterilisation program types .....                        | 56        |
| 4.3.2 Cycle B Diagram .....                                    | 56        |
| 4.3.3 Cycle S Diagram.....                                     | 56        |
| 4.4 Sterilisation program specifications .....                 | 57        |
| 4.5 Putting the material to be sterilised in the chamber ..... | 57        |
| 4.5.1 Preparing the material.....                              | 57        |
| 4.6 Filling and draining water.....                            | 58        |
| 4.6.1 Filling the feed tank.....                               | 58        |
| 4.6.2 Draining the discharge tank .....                        | 58        |
| 4.7 Test programs.....                                         | 59        |
| 4.7.1 Test program types.....                                  | 60        |
| 4.7.2 Vacuum test (leakage test).....                          | 60        |
| 4.7.3 Bowie & dick test.....                                   | 60        |
| 4.7.4 Helix test .....                                         | 60        |
| 4.7.5 Test diagrams .....                                      | 60        |

|                                      |    |
|--------------------------------------|----|
| 4.8 Settings .....                   | 61 |
| 4.8.1 User settings .....            | 61 |
| 4.8.2 Date and time settings.....    | 61 |
| 4.8.3 Printer settings .....         | 62 |
| 4.8.4 Language settings .....        | 62 |
| 4.8.5 Measurement unit settings..... | 62 |
| 4.8.6 Water settings .....           | 62 |
| 4.8.7 Connections .....              | 63 |
| 4.8.8 Ethernet connection .....      | 63 |
| 4.9 Data output .....                | 66 |

## 5 Routine maintenance .....

|                                                                        |    |
|------------------------------------------------------------------------|----|
| 5.1 General table .....                                                | 69 |
| 5.1.1 Cleaning the sterilisation chamber and internal components ..... | 69 |
| 5.1.2 Cleaning the door cover .....                                    | 69 |
| 5.1.3 Cleaning the external front panel.....                           | 69 |
| 5.1.4 Cleaning the door seal.....                                      | 69 |
| 5.1.5 Cleaning the filter in the chamber.....                          | 69 |
| 5.1.6 Replacing the Bacterial or Bio X filter .....                    | 69 |
| 5.1.7 Replacing the door seal.....                                     | 69 |
| 5.1.8 Cleaning the dust filter (if provided) .....                     | 69 |
| 5.1.9 Lubricating and cleaning the screw and nut screw .....           | 70 |
| 5.1.10 Clean the feed tank .....                                       | 70 |
| 5.2 Maintenance.....                                                   | 70 |
| 5.3 Periodic process effectiveness and safety checks .....             | 70 |

## Troubleshooting

|                                                             |           |
|-------------------------------------------------------------|-----------|
| <b>6 Troubleshooting guide.....</b>                         | <b>71</b> |
| 6.1 Error Codes .....                                       | 71        |
| 6.2 Alarm Codes.....                                        | 72        |
| 6.3 Messagges.....                                          | 72        |
| 6.4 Reset procedures .....                                  | 73        |
| 6.5 Troubleshooting .....                                   | 73        |
| 6.6 Message solutions .....                                 | 76        |
| <b>7 Technical specifications .....</b>                     | <b>76</b> |
| 7.1 Water specifications (distilled or demineralized) ..... | 76        |
| <b>8 Spare parts codes .....</b>                            | <b>76</b> |

## ! Important information

### SYMBOLS

#### HAZARD

Paragraphs marked with this symbol contain instructions that you must follow carefully to avoid damaging the device, injuring the operator and possibly the patient.

#### WARNING

These instructions warn that you must take particular care to avoid situations that could damage the device.

#### PROHIBITION

This icon highlights what you must not do, to avoid damaging the device.

#### TIPS

With this icon, you are given information to allow you to use the device more effectively.

#### WARNING

General warning symbol.

#### WARNING

Risk or hazard symbol. Consult the manual.

### SYMBOLS SHOWN ON THE PRODUCT

0051

European conformity marking (CE) with identification number of the Notified Body for Class IIa, IIb, III medical devices. Notified Body No. 0051: IMQ S.p.A., Italy

OFF symbol. Located on the mains switch (Fig. 1 §2).

Indicates the manufacturer serial number with which a specific medical device can be identified.

Consult the user manual or the service manual.

After its useful life, the device must be collected separately from other waste according to the differentiated waste collection regulations. Illegal disposal of the product will result in application of fines defined by the individual legislations.

When the device is in operation, there must always be a tray holder and a tray positioned on the bottom shelf.

This symbol is located in proximity of the water loading point to be used in the event of an emergency.

This symbol is located in proximity of the water draining point.

This symbol is located in proximity of the water loading point.

This symbol is located in proximity of the waste water draining point.

This symbol is located in proximity of the points on the machine where there is a potential risk of burns.

This symbol is located next to the USB connection point to the device.

Indicates the medical device manufacturer as **defined in the EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC**. This symbol is accompanied by the name and address of the manufacturer.

Indicates the catalogue number of the manufacturer with

### PACKAGING SYMBOLS


HIGH

DO NOT WET

DO NOT USE HOOKS

STORAGE TEMPERATURE  
TRANSPORT CONDITION

P

500 hPa ÷ 1060

ATMOSPHERIC PRESSURE

FRAGILE

DO NOT ROLL

MAX SUPPORTED  
WEIGHT

RH

10% ÷ 90%

RELATIVE HUMIDITY

RECYCLABLE CARDBOARD

### PACKAGING

We recommend keeping the packaging for handling the device in the future. In the event of disposal, separate the various parts in accordance with the differential collection regulations.

# 1 Safety standards

Follow all instructions in this manual.

The autoclave may only be used by personnel specially trained for the sterilization process in medical practices, podiatry practices and beauty centers under the responsibility of the medical director and the supervisor for the sterilization process.

The autoclave may only be used by personnel who have been instructed on the contents of these instructions for use. The medical facility in which the autoclave is installed is responsible for training in the operation and maintenance of the product. It must keep training records and check that the personnel have understood the training. Assembly, upgrades, maintenance and repair may only be carried out by Dürr Dental or by an agency authorized by Dürr Dental for this purpose. The operator is responsible for ordering service only from authorized agencies. There are no procedures by which the user can check the correct function of the unit himself.

## 1.1 Intended users

- medical staff (dentists, chiropodists)
- dental hygienists
- employees in beauty salons

### USER PROFESSIONAL QUALIFICATIONS

Degree in medicine

Dental hygienist

### MINIMUM SKILLS

Those required by the professional qualification

Course on use based on this manual for unqualified users

### LANGUAGE COMPREHENSION

Those acquired for the professional qualification

### EXPERIENCE

That required to carry out the profession

### POSSIBLE USER HANDICAPS

Full use of an upper limb is required

Eyesight compatible with the profession

## 1.2 Intended use

It is a small steam autoclave for sterilising tools, rotary instruments, reusable medical devices and porous materials (e.g. fabrics) using saturated steam at a temperature of 121°C or 134°C.

## 1.3 General safety warnings

The product must be used for its intended use in accordance with the provisions of this manual.

Any use that does not comply with the instructions in the instruction manual could compromise the safety level of the device.

Do not make any modification to the equipment without the written authorization by the manufacturer.

The user is responsible for installation, training users, maintenance and any legal requirements related to the device, including the periodic checks described in this manual.

DÜRR DENTAL shall not be liable for harm to things and/or people caused by failure to use the product correctly as described in this manual. Failure to comply with the provisions of this manual will void any liability of the manufacturer and the product warranty.

The device discharges air from a small hole at its bottom rear. Do not obstruct this vent (**17 fig. 2**).

Do not place the autoclave on an unstable table or surface.

### Electrical hazard warnings

The device must be connected to the electrical mains. The electrical system must comply with the IEC 364-1 standard and the "National regulations for installation of electrical systems in premises intended for medical purposes".

Check that the mains supply matches the power supply voltage indicated on the plate on the autoclave rear panel. If in double, consult your dealer.

The plug on the autoclave has an earth connection. Make sure that you use an earthed socket and position the device where it is not difficult to disconnect the plug.

If you use an extension cable, make sure that the cable type is adequate for the equipment consumption.

Do not attempt to carry out maintenance on the autoclave when it is connected to the power. Disconnect the power cable from the mains before working on the machine. Maintenance or repair operations must be carried out by persons that have obtained a training certificate issued by the manufacturer or Dürr Dental in the capacity of Exclusive Distributor.

**Do not carry out maintenance work other than that described in this manual. Any work not indicated in this manual could compromise equipment safety. Contact DÜRR DENTAL technical service for any work that is not specified.**

### Equipment explosion hazard warnings

Do not use any acids or corrosive substances to clean any parts of the autoclave. Specifically, do not use hydrochloric acid or any detergent that contains chlorine. Using these substances may compromise the mechanical integrity of the equipment and create an explosion hazard.

The pressure chamber is protected by a safety valve: annually check that it is working correctly.

**Do not carry out maintenance work other than that described in this manual. Any work not indicated in this manual could compromise equipment safety. Contact DÜRR DENTAL technical service for any work that is not specified.**

### Contamination hazard warnings

The autoclave must be installed in a suitable hygienically controlled environment.

**Clean the chamber before using the device for the first time. The autoclave works with distilled or demineralised water (see technical specifications in chapter 7.1 or described in EN 13060 Appendix C). The device would not work properly with any other substance. Consult DÜRR DENTAL technical service if an unsuitable liquid is put inside.**

Do not use toxic substances during the sterilisation process. Do not touch the contents of the discharge tank.

Do not put substances or products that contain heavy metals (lead, mercury, bromine, hexavalent chromium) into the chamber. If released, these products could contaminate the autoclave plumbing circuit and render it unusable.

### NEVER REUSE WATER

**Do not carry out maintenance work other than that described in the manual.**

Only put heat-resistant products that do not release hazardous substances into the sterilisation chamber. Always consult the manufacturer's technical specifications to ensure that the products can be sterilised.

All items must be decontaminated and thoroughly cleaned and dried before being sterilised.

We recommend using chemical indicators to confirm the sterilisation process.

### Fire or overheating hazard warnings

Only put products that are resistant to steam sterilisation inside, in accordance with the manufacturer's provisions. Do not place the device on a flammable surface.

The device is not suitable for use in oxygen-rich or potentially explosive atmospheres.

Do not block the air vents on the equipment casing during installation or operation.

Comply with the minimum distance requirements in detail if the autoclave is recessed.

#### Operation failure hazard warning

If an error message appears while a sterilisation cycle is running, it will be necessary to repeat the cycle. Accidentally dropping the device may deform it and compromise proper operation. If this happens, have a full check done by technical support. If there are instruments that are joined together, they must be divided or located in the most ventilated and spacious place possible. Lab coats or other reusable fabrics must be washed and dried after use and before sterilisation to remove any organic material and increase the "life" of the fabric by restoring its natural water content (i.e. humidity level). Comply with the load limits defined in this manual for each type of sterilisation load.

#### Crushing of body parts hazard warnings

The equipment has motorised door closure. Only hold the door by the handle, or push it from the outside. Be careful of your finger/hand position while closing the door. Do not leave fingers or other body parts between the door and the machine while it is closing.

#### Burn hazard warnings

To prevent possible burns, extract the trays with the dedicated pincers provided. Before cleaning any surface, make sure that it has cooled down.

#### Infect hazard warnings

The autoclave performs a sterilisation process. However, the user must guarantee the effectiveness of the sterilisation process by following the guidelines as well as the data provided by the autoclave. Therefore, strictly follow the protocols defined in this manual to ensure that the process continues to be effective.

#### Warnings against the risk of operator infection

Despite the water in the drain tank undergoing the sterilization process, there is still a risk of biological contamination. It is therefore prohibited to touch the drain tank water. For disinfection follow the recommended procedures in Section 5.

#### Note

The device, when on, must always have the tray holder and tray on the lower shelf.

### 1.4 Operating and transport environmental conditions

In its packaging, the equipment can withstand the following environmental conditions for no more than 15 weeks:

Ambient temperature from -20°C to +70°C

Relative humidity from 10% to 90%

Atmospheric pressure from 500 to 1060 mBar.

The appliance must be used in the following environmental conditions:

- internal use
- altitude up to 3000 m
- Atmospheric pressure from 700 hPa to 1060 hPa
- temperature from 5°C to 40°C
- max relative humidity 85%
- max mains voltage variation  $\pm 10\%$
- installation category (overvoltage category) II
- pollution degree rating 2
- ambient lighting 500 lx

### 1.5 Environmental disposal instructions


Pursuant to directives 2012/19 EC and 2011/65 EC on waste electrical and electronic equipment and the restriction of the use of certain hazardous substances in electrical and electronic equipment, waste cannot be disposed of as urban waste, but must be collected separately.

When purchasing new equipment of an equivalent type on a one-for-one basis, the equipment that has reached the end of its life must be returned to the dealer for disposal. Regarding reuse, recycling or the other forms of waste recovery mentioned above, the manufacturer must carry out the steps defined by the individual national legislation.

Adequate separate collection for future use of the equipment assigned for recycling, treatment and environmentally friendly disposal helps to prevent possible negative effects on the environment and on human health, and promotes recycling and/or reuse of the materials that make up the equipment.

The crossed-out bin symbol (see symbols on page 38) on the device or on its packaging indicates that the product must be disposed of separately from other waste at the end of its life.

## 2 Equipment description


**Fig. 1**

1. Touch screen display
2. Bio-x filter
3. Motorised closure screw
4. USB port
5. Main switch
6. Feed tank filling connection
7. Discharge tank draining connection
8. Feed tank draining connection
9. Manual feed tank emergency filling connection
10. Nut screw for motorised closure
11. Door cover
19. Dust filter
21. Tray holder symbol
22. Loading tank opening
23. Drain tank opening


**Fig. 2 – Rear view**

12. Power supply
13. RS-232 port for external printer
14. LAN port
15. Discharge tank drain connection for main drain system
16. Discharge tank vent
17. Feed tank vent
18. Solenoid valve for filling the feed tank from the mains water supply (Allowable pressure from 20 kPa to 250 kPa)

**Fig. 2 and Fig. 3**

20. Data plate

## 2.1 Variants

The following autoclave variants are available:

| SERIES     | HC 50                           |                   |
|------------|---------------------------------|-------------------|
| TYPE       | Hygoclave 50                    | Hygoclave 50 Plus |
| REF        | 604 910 XXXX                    | 604 920 XXXX      |
| Values XXX |                                 |                   |
| 0001       | Only type B, 230V, 50/60 Hz     |                   |
| 0002       | Only type B, 240V, 50 Hz        |                   |
| 0003       | CH, Only type B, 230V, 50/60 Hz |                   |
| 0004       | Only type B, 220V, 50/60 Hz     |                   |
| 0005       | Only type S, 230V, 50/60 Hz     |                   |
| 0006       | Only type S, 240V, 50 Hz        |                   |
| 0007       | CH, Only type S, 230V, 50/60 Hz |                   |
| 0008       | Only type S, 220V, 50/60 Hz     |                   |
| 0009       | Type B+S, 230V, 50/60 Hz        |                   |
| 0010       | Type B+S, 240V, 50 Hz           |                   |
| 0011       | CH, Type B+S, 230V, 50/60 Hz    |                   |
| 0012       | Type B+S, 220V, 50/60 Hz        |                   |
| 0013       | UK, Only type B, 230V, 50/60 Hz |                   |
| 0014       | UK, Only type S, 230V, 50/60 Hz |                   |
| 0015       | UK, Type B+S, 230V, 50/60 Hz    |                   |
| 0025       | AU, Type B+S, 230 V, 50/60 Hz   |                   |

## 2.2 Main features

The vacuum pump ensures an optimal sterilisation process in which the steam interacts properly with the load being processed, including any cavities.

### 2.2.1 Sterilisation programs

- **Type B**
  - 134°C **Standard**
  - 134°C **Standard +**
  - 134°C **Prion**
  - 134°C **B-Flash**
  - 121°C **Standard**
- **Special cycle**
  - Drying only
- **Type S**
  - 134°C **Universal S**
  - 121°C **Softprogramm**

### 2.2.2 Test programs

- **Vacuum test:** checks the ability to hold a vacuum
- **Bowie & Dick test:** checks steam penetration into porous loads
- **121 or 134°C Helix test:** checks steam penetration into hollow loads

The autoclave is available with the following storage devices:

- USB port (USB peripheral not supplied).
- RS-232 port (for DÜRR DENTAL printer, NOT supplied).

### 2.2.3 Water filling systems

**Automatic:** with the aid of a pump inside the autoclave.

**Mains connection** possibility of connecting a demineralisation system, to increase autonomy. It must be installed by an authorised technician.

**Manual:** as an alternative, the device is equipped with a manual emergency filling system at the FRONT (9 fig. 1).

### 2.2.4 Automatic sterilisation parameter control

During sterilisation, the pressure, temperature and time are continuously monitored by an automatic process control system.

### 2.2.5 Motorised closure

The autoclave has a motorised closure system to ensure a good seal. It is based on a screw and nut screw. Closure works with a pre-locking mode that activates when the door is shut. The door is only closed fully when the cycle starts.

### 2.2.6 Bio-X filter

This injects sterile air into the chamber during the drying and pressure balancing phases.

### 2.2.7 USB port

This port can be used by medical staff for data download. Exclusively use a USB Key from version 1.1 with a capacity of less than 8 GB.

### 2.2.8 Serial port - RS232

Port dedicated exclusively to connection of the Hygoprinter Easy printer. Use the dedicated serial cable provided with the printer. Use only the dedicated Hygoprint Easy printer (order number 6048100167).

### 2.2.9 Touch screen display

This manages the autoclave fully through a menu system.

### 2.2.10 LAN connection

This port is only for connecting a LAN. Protocol used: 100Base-T. Connect only an Ethernet cable.

### 2.2.11 Safety devices

The steriliser is equipped with the following safety devices:

#### - Mains fuses

These protect the equipment from faults in the heating elements by disconnecting the electric power supply.

#### - Electronic circuit protection fuses

These prevent faults in the transformer primary circuit and low-voltage loads by disconnecting one or more low-voltage circuits.

#### - Safety valve

This prevents any overpressure in the sterilisation chamber by releasing steam to restore a safe pressure. It opens at a pressure of 350kPa ( $\pm 10\%$ ).

#### - Steam generator thermal switch

This stops the steam generator from overheating by disconnecting its electrical power supply. It trips at a temperature of  $220^{\circ}\text{C} \pm 7^{\circ}\text{C}$ .

#### - Chamber thermal switch

This protects the pressure vessel heating element from overheating by disconnecting its electric power supply. It trips at a temperature of  $150^{\circ}\text{C} \pm 4^{\circ}\text{C}$ .

#### - Door lock mechanism safety microswitch

This provides feedback when the door and the door lock system are not in the correct position or when the door lock mechanism is not working correctly.

#### - Pressure switch

This prevents the door from being opened accidentally during the program when the machine is pressurised (even in the event of a blackout).

#### - Overflow system

This prevents water from flowing out onto live parts if the level sen-

sor does not work by draining excess water to the outside.

#### - Built-in sterilisation process assessment system

This uses a microprocessor to continuously monitor the sterilisation process parameters. In the event of a fault, it immediately stops the program and generates alarms.

#### - Steriliser operation monitoring

When the machine is powered, this monitors all of its important parameters in real time. In the event of a fault, it generates alarm messages and, if necessary, stops the cycle.

### 2.3 Accessories included in the medical device certification

Only the following accessories are included in the medical device certification:

- The tray support
- 5 trays
- 1 tray extraction pincer
- 1 filling hose
- 1 drain hose
- 1 funnel
- 1 bottle of lubricating oil

Any other accessory mentioned in this manual does not form part of the CE 0051 certification.

#### TYPE TESTS CARRIED OUT ON AUTOCLAVES

Type test carried out in accordance with the EN13060 standard

Type B cycles

|                                             | Cycles B | Cycles S |
|---------------------------------------------|----------|----------|
| Pressure dynamic in sterilisation chamber   | X        | X        |
| Air leakage                                 | X        | X        |
| Empty chamber                               | X        | X        |
| Solid load - Double wrapping                | X        | X        |
| Small porous load<br>Double wrapping        | X        | X        |
| Full porous load<br>Double wrapping         | X        | -        |
| Narrow lumen element                        | X        | X        |
| Drying, solid load<br>Double wrapping       | X        | X        |
| Drying, full porous load<br>Double wrapping | X        | -        |

## 3 Installation

**THE DEVICE MUST BE INSTALLED BY QUALIFIED TECHNICIANS.**

- Check that the installation spaces comply with the indications given in the dimensions. To recess the autoclave inside a cabinet, there must be enough space all around (20 mm sideways, 50 mm at the top, 30 mm at the back) it to ensure effective ventilation, and a wide opening at the rear for the power cable. Make sure that there is sufficient air flow to properly cool the heat exchanger.

- Install the steriliser on a perfectly horizontal flat surface and make sure that it is strong enough to support the weight of the device (approximately 60 kg); To prevent contact with water or other liquids that could cause short circuits and/or potential danger to the operator, do not install the steriliser too close to basins, sinks or similar places;

- Do not install the steriliser in excessively damp or poorly ventilated places;


- Do not install the machine in environments where there is gas, flammable vapours and/or explosives;

- Position the equipment so that the mains switch is easily accessible.

- Install the equipment so that the power cable is not bent or pinched. It must be free to move to the power socket.

- Install the equipment so that any external filling/drain hoses are not bent or pinched. They must be free to move to the drain tank.

### External dimensions

**There are two lifting points on the autoclave base.**

**When handling the equipment briefly, it must be done by two people; the lifting points (for straps) are located on the equipment feet. For longer handling, insert the pallet provided under the equipment.**

Free the autoclave from the packaging.

- Position the autoclave on the support (table or cabinet) using straps.

- Make sure that the supporting surface is not made of flammable material.

- Remove the straps and keep them for future handling.

- **Use a spirit level to check that the support base is level in both directions, otherwise level the autoclave using its adjustable front feet.**

- **Check that the power socket used is earthed.**

- **Check that the mains voltage matches that on the autoclave plate.**

- Insert the autoclave plug, making sure that the cable is not pinched and that it can be unplugged easily.

- The plug is the mains means of connection to the mains power, so it must be easy for the operator to access and use it.

### Electrical connections

In accordance with the laws and/or regulations in force, the steriliser must be connected to an earthed power socket with adequate current capacity for the equipment consumption (see the data plate).

The socket must be properly protected by a residual current circuit breaker with the following specifications:

- Nominal current  $I_n$  16 A

- Differential current  $I_{\Delta n}$  0.03 A

**CAUTION:** The manufacturer shall not be liable for damage caused by installing the steriliser with inadequate and/or unearthed electrical systems.

**NOTE:** Always connect the power cable directly to the power socket.

Do not use extension cables, adapters or other accessories.

### 3.1 FIRST POWER-UP

- Press the main switch (5 Fig. 1). The green light on the switch will turn on to indicate that the device is powered.
  - 1 - Wait for the main screen to load
  - 2 - Select the desired language.
  - 3 - "Warranty Activation" screen.
  - 4 - Home page loading following the choices made during warranty activation.
- Select the "07 open door" icon to open the door and extract the supplied accessories from the chamber.

#### EXAMPLE OF LANGUAGES SCREEN

#### HOME PAGE SCREEN

- 01 "Sterilisation cycles" menu
- 02 "Test" menu
- 03 "Settings" menu
- 04 "Data Output" menu
- 05 Last sterilisation cycle run
- 06 Last test run
- 07 Unlock door
- 08 Water filling

**"WARRANTY ACTIVATION" SCREEN** on the website [www.DÜRR DENTAL.it](http://www.DÜRR DENTAL.it):  
The screen will appear every 24 hours for 30 days. Please register your device on Dürr Dental's website ([www.duerrdental.com](http://www.duerrdental.com)) to benefit from an additional year of warranty free of charge. When you press the "A" button you will see the "Registration Confirmation" screen.


#### Check that the accessories include:

- Documentation (instruction manual)
- The tray support
- 5 trays
- 1 tray extraction pincer
- 1 filling hose
- 1 drain hose
- 1 funnel

#### WARNING

- Insert the tray holder as shown in the figure
- The tray holder must be inserted into the chamber with foot "A" facing down and positioned at the bottom of the chamber. The position of foot "A" must be respected both when using 3 trays and 5 trays
- The tray holder must always have rubber protection "B" mounted on the eight supports

If the fresh water tank is filled neither manually nor semi-automatically (using an external canister), the autoclave **MUST NOT** be connected directly to the mains water supply. In this case (automatic filling), the fresh water supply must be provided by an intermediate water treatment system, which must be able to provide water according to the required specifications (see 7.1 or EN 13060 Appendix C). The water treatment system must prevent water backflow to the mains.

#### Validate the installation in accordance with the following table:

- Vacuum test
- Autoclave automatic test
- Autoclave instrumentation calibration check
- Overtemperature trip test
- Thermometric test for small loads (\*)
- Drying test for small loads (\*)
- Steam penetration test for porous loads (\*)
- Thermometric test for solid loads (\*)
- Drying test for solid loads (\*)
- Steam penetration test for solid loads (\*)
- Thermometric test for a load chosen by the user
- Microbiological test for a load chosen by the user

(\*) To be carried out with cycles of 121°C and 134°C

## 4 Operating instructions

Operate the autoclave as follows:


- Press the main switch (**5 fig. 1**). The green light on the switch will turn on to indicate that the device is powered.
- Wait for the main screen to load

### USING THE MENUS

- 01 Sterilisation programs
- 02 Test programs
- 03 Settings
- 04 Data output
- 05 Last cycle run
- 06 Last test run
- 07 Open door
- 08 Water filling

## 4.1 SYMBOLS

| Description             | Icon                                                                                | What it does/what it means                                                                                                           |
|-------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| Screen name             |    | Indication icon only. Cannot be selected. Shows the name of the selected screen.                                                     |
| Time                    |    | Indication icon only. Cannot be selected. Shows the name of the selected screen.                                                     |
| Date                    |    | Indication icon only. Cannot be selected. Shows the name of the selected screen.                                                     |
| Back                    |    | Returns to the previous screen.                                                                                                      |
| Home                    |    | Returns to the home screen.                                                                                                          |
| Confirm                 |    | Confirms certain kinds of actions.                                                                                                   |
| Cycles                  |    | Goes to the cycles screen where you can select the desired cycle.                                                                    |
| Test                    |    | Goes to the tests screen where you can select the desired test.                                                                      |
| Settings                |    | Goes to the machine settings screen.                                                                                                 |
| Data output             |    | Goes to the data output screen.                                                                                                      |
| Open door               |   | Opens the door.                                                                                                                      |
| Water filling           |  | Goes to the water filling screen.                                                                                                    |
| 134°C Standard cycle    |  | Goes to the 134°C standard cycle screen, which displays all of the cycle specifications.                                             |
| 134°C Standard+ cycle   |  | Refers to the Cycle 134°C Standard+ screen (with 4 instead of 3 fractionations levels), in which all cycle properties are specified. |
| 121°C Standard cycle    |  | Goes to the 121°C standard cycle screen, which displays all of the cycle specifications.                                             |
| 134°C Prion cycle       |  | Goes to the 134°C Prion cycle screen, which displays all of the cycle specifications.                                                |
| 134°C B Flash cycle     |  | Goes to the 134°C B Flash cycle screen, which displays all of the cycle specifications.                                              |
| Drying cycle            |  | Goes to the drying cycle screen. This is considered a special cycle to be used when an additional drying cycle is needed.            |
| Programmed start        |  | Programs the cycle start time. The set time is displayed beside the corresponding icon.                                              |
| Increment drying        |  | Increments the drying time. The value entered will update the drying time in the cycle.                                              |
| Next                    |  | Increases the value in the relative field by one or goes to the next option.                                                         |
| Previous                |  | Decreases the value in the relative field by one or goes to the previous option.                                                     |
| Vacuum Test             |  | Goes to the Vacuum Test screen, which displays all of the cycle specifications.                                                      |
| 134°C Helix Test        |  | Goes to the 134°C Helix Test screen, which displays all of the cycle specifications.                                                 |
| 134°C Bowie & Dick Test |  | Goes to the 134°C Bowie & Dick Test screen, which displays all of the cycle specifications.                                          |

| Description              | Icon                                                                                | What it does/what it means                                                                                                                                                                                                                                 |
|--------------------------|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 121°C Helix Test         |    | Goes to the 121°C Helix Test screen, which displays all of the cycle specifications.                                                                                                                                                                       |
| Programmed start         |    | Programs the test start time.                                                                                                                                                                                                                              |
| Next                     |    | Increases the value in the relative field by one unit or goes to the next option.                                                                                                                                                                          |
| Prevoius                 |    | Decreases the value in the relative field by one unit or goes to the previous option.                                                                                                                                                                      |
| User settings            |    | Goes to the user setting screen.<br>You can choose from: add user, delete user, edit/change user.<br>This section can only be accessed by the administrator user "Admin", who will be asked to enter a PIN.                                                |
| Date/time settings       |    | Goes to the data/time settings screen. Sets the time (hours, minutes) and date (dd/mm/yyyy).                                                                                                                                                               |
| Language settings        |    | Goes to the language settings screen. Selecting a different language automatically returns to the setting screen with all the fields updated in the selected language.                                                                                     |
| Mesurement unit settings |    | Goes to the measurement unit settings screen.<br>Changes the pressure measurement unit.                                                                                                                                                                    |
| Connectivity settings    |    | Goes to the connectivity settings screen.                                                                                                                                                                                                                  |
| Water parameter settings |   | Goes to the water parameter settings screen.<br>Changes the water filling type (manual at the front or water mains at the rear).<br>Displays the water quality. There are three levels: green, yellow and red, with the cursor on the corresponding level. |
| Printer settings         |  | Goes to the pinter settings screen.<br>Sets manual or automatic printing.                                                                                                                                                                                  |
| Service                  |  | Goes to the service action, after entering a PIN.<br>The service section has dedicated graphics to show users where they are.                                                                                                                              |
| Next                     |  | Increases the value of the relevant field by one unit or move to the next option.                                                                                                                                                                          |
| Prevoius                 |  | Decreases the value of the relevant field by one unit or move to the previous option.                                                                                                                                                                      |
| ENGLISH                  |  | Sets the language for all menus to ENGLISH.                                                                                                                                                                                                                |
| ITALIAN                  |  | Sets the language for all menus to ITALIAN.                                                                                                                                                                                                                |
| FRENCH                   |  | Sets the language for all menus to FRENCH.                                                                                                                                                                                                                 |
| GERMAN                   |  | Sets the language for all menus to GERMAN.                                                                                                                                                                                                                 |
| SPANISH                  |  | Sets the language for all menus to SPANISH.                                                                                                                                                                                                                |
| TURKISH                  |  | Sets the language for all menus to TURKISH.                                                                                                                                                                                                                |
| POLISH                   |  | Sets the language for all menus to POLISH.                                                                                                                                                                                                                 |
| PORTUGUESE               |  | Sets the language for all menus to PORTUGUESE.                                                                                                                                                                                                             |
| CZECH                    |  | Sets the language for all menus to CZECH.                                                                                                                                                                                                                  |
| HUNGARIAN                |  | Sets the language for all menus to HUNGARIAN.                                                                                                                                                                                                              |

| Description                       | Icon                                                                                | What it does/what it means                                                                                                                            |
|-----------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| Croatian                          |    | Sets the language for all menus to CROATIAN.                                                                                                          |
| Greek                             |    | Sets the language for all menus to GREEK.                                                                                                             |
| Ethernet settings                 |    | Goes to the Ethernet settings screen, where you can set the network parameters (IP address, subnet mask, gateway).                                    |
| Add users                         |    | Goes to the user data entry screen (name, surname, PIN).                                                                                              |
| Edit users                        |    | After you have selected which user to edit, it goes to the user edit scan (name, surname, PIN).                                                       |
| Delete users                      |    | After you have selected which user to delete, it asks for confirmation.                                                                               |
| Admin user                        |    | Goes to the Admin data (name, surname, PIN).                                                                                                          |
| Generic user                      |    | From the edit section, it goes to the user data (name, surname, PIN).<br>From the deleted section, it determines which user is to be deleted.         |
| Alphanumeric keypad               |    | Allows you to enter letters, numbers and symbols .<br>Eg. To select the letter C, press the corresponding key 4 times (2, A, B, C).                   |
| Usb                               |    | Goes to the screen whre you can load cycles onto a USB pen drive. After inserting the USB drive, you can select the cycles that you want to transfer. |
| Print                             |    | Goes to the screen whre you can print reports, labels and barcodes.<br>The icon will be disabled if the printer is not connected.                     |
| Info                              |  | Goes to the screen with the machine manufacturing information: SN, PN, installation date, configuration, FW, no. of cycles, website.                  |
| Barcode                           |  | Selects barcode printing. Once selected, a menu appears from which you can choose the number of labels to print. Printing starts when you confirm.    |
| Report                            |  | Selects and activates cycle report printing.                                                                                                          |
| Label                             |  | Selects label printing. Once selected, a menu appears from which you can choose the number of labels to print. Printing starts when you confirm.      |
| Increase                          |  | Increases the value in the relative field by ten.                                                                                                     |
| Decrease                          |  | Decreases the value in the relative field by ten.                                                                                                     |
| STOP water filling                |  | Stops the muanual water filling pump.                                                                                                                 |
| START water filling               |  | Starts the muanual water filling pump.<br>The button is disabled when filling from the water mains.                                                   |
| Previous                          |  | Goes to the previous option in the scrolling menu.                                                                                                    |
| Next                              |  | Goes to the next option in the scrolling menu.                                                                                                        |
| Next                              |  | Increases the value in the relative field by one or goes to the next option.                                                                          |
| Previous                          |  | Dencreases the value in the relative field by one or goes to the previous option.                                                                     |
| Programmed start icon             |  | Indication icon only. Cannot be selected.<br>Shows the start time of the selected cycle.                                                              |
| Pressure indicator                |  | Indication icon only. Cannot be selected.<br>Shows the pressure in the selected cycle.                                                                |
| Maximum instrument load indicator |  | Indication icon only. Cannot be selected.<br>Shows the maximum instrument load for the selected cycle.                                                |


| Description                    | Icon                                                                                | What it does/what it means                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|--------------------------------|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Remaining time indicator       |    | Indication icon only. Cannot be selected.<br>Shows the remaining drying time in the cycle (to be checked) .                                                                                                                                                                                                                                                                                                                                                                                   |
| Temperature indicator          |    | Indication icon only. Cannot be selected.<br>Shows the temperature in the selected cycle.                                                                                                                                                                                                                                                                                                                                                                                                     |
| Maximum textile load indicator |    | Indication icon only. Cannot be selected.<br>Shows the maximum textile load for the selected cycle.                                                                                                                                                                                                                                                                                                                                                                                           |
| Door lock indicator            |    | Indication icon only. Cannot be selected.<br>Shows that the door is closed and locked.                                                                                                                                                                                                                                                                                                                                                                                                        |
| Cycle count indicator          |    | Indication icon only. Cannot be selected.<br>Shows the cycle counter.                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Stop cycle/test                |    | Stops a cycle/test at any time.<br>Press the icon for at least 3 seconds to activate it. After 3 seconds, the cycle/test will stop and the machine is put in a safe state.                                                                                                                                                                                                                                                                                                                    |
| Cycle completed confirmation   |    | Indication icon only. Cannot be selected.<br>Appears when the cycle has ended correctly.                                                                                                                                                                                                                                                                                                                                                                                                      |
| Cycle error                    |    | Indication icon only. Cannot be selected.<br>Appears when the cycle has been stopped by an error or manually by the operator.                                                                                                                                                                                                                                                                                                                                                                 |
| Cycle progress bar             |    | The cycle progress bar advances according to the following phases:<br><ol style="list-style-type: none"> <li>1. Preheating phase</li> <li>2. Vacuum phase 1</li> <li>3. Pressure ramp-up phase 1</li> <li>4. Vacuum phase 2</li> <li>5. Pressure ramp-up phase 2</li> <li>6. Vacuum phase 3</li> <li>7. Pre-sterilisation phase</li> <li>8. Balancing phase</li> <li>9. Sterilisation phase</li> <li>10. Discharge phase</li> <li>11. Drying phase</li> <li>12. Pressure balancing</li> </ol> |
| Sterilisation complete         |  | Indication icon only. Cannot be selected. Appears when the sterilisation phase is complete (not the cycle). Allows the user to stop the cycle to remove the instruments, skipping the drying phase.                                                                                                                                                                                                                                                                                           |
| Feed tank maximum level        |  | Appears when the water in the feed tank reaches the maximum level. Water filling stops automatically.                                                                                                                                                                                                                                                                                                                                                                                         |
| Feed tank intermediate level   |  | Appears when the water in the feed tank is between the minimum and maximum levels.                                                                                                                                                                                                                                                                                                                                                                                                            |
| Feed tank minimum level        |  | Appears when the water in the feed tank is close to the minimum level.                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Water quality                  |  | Indication icon only. Cannot be selected.<br>Shows the water quality bar.                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Water quality bar              |  | The cursor position indicates the water quality.<br>Green: quality OK<br>Yellow: critical quality area → message generated<br>Red: quality not OK → message generated and machine stops.                                                                                                                                                                                                                                                                                                      |

## 4.2 Menu architecture

- Home
  - Sterilisation programs
  - Last sterilisation program run
  - Test programs
  - Last test program run
- Settings
  - User management
  - Date and time
  - Printer
  - Language
  - Measurement unit
  - Water parameters
  - Connectivity
  - Service
- Data output
  - Usb
  - Printer
  - Info
- Open door
- Water filling

### 4.3 Selecting sterilisation programs

- 01 Select the "Sterilisation cycles" menu
- 02 Alternatively, directly select "Last cycle run" and go to step 04
- 03 Select the desired cycle
- 04 Screen with all the parameters of the selected cycle
- 05 (Optional) You can open them door to insert other instruments
- 06 (Optional) Select the programmed start button to set the start time
- 07 (Optional) Select the increment drying button to increase the cycle drying time
- 08 Start the cycle
- 09 (If active) enter the user PIN.
- On confirming, the cycle will start

- 10 You can use the green bar to monitor the cycle progress
- 11 When the cycle finishes, press the confirm button to return to the Home screen. If the Cycle Validation Protocol Management is enabled, at the end of the cycle, when you confirm the selection by, the system shows the Cycle Validation Protocol screen.
- 12/13 You can stop the cycle manually by holding done the "Stop cycle" buttonh for 3 seconds. Wait for the machine to be put in a safe condition and press the confirm button (13) for 2 seconds and then press Home to return to the home screen.
- 14/15 Cycle Validation Protocol: The selection of the Batch, Load and/or Indicator validation options generates the relative information in the cycle report: "Batch Ok" - "Batch Fail" "Indicator OK" - "Indicator Fail" - "No Indicator". "Load OK" - "Load Fail" - "No Load" \* The confirmation button (15) will be "Operational" when the user makes a selection.
- 16 For the validation of Batch, Load and/or Indicator, the user must enter his PIN. When the confirmation button is pressed (16), the system generates the Cycle Report PDF with the identification of the User who performed the cycle validation.

### 4.3.1 Sterilisation program types

When the machine is cold, the preheating time is approximately 10 minutes. Selecting the “Programs” submenu accesses the list of sterilisation cycles that the autoclave can run:

#### B CYCLES

- **134°C Standard cycle:** 5.5 min. sterilization plus drying phase.
- **134°C Standard + cycle:** 5.5 min. sterilization plus drying phase. 4 vacuum
- **134°C Prion cycle:** 20 min. sterilization plus drying phase. 4 vacuum phases
- **134°C B-Flash cycle:** 5.5 min. sterilization plus drying phase.
- **121°C Standard cycle:** 20.5 min. sterilization plus drying phase.

#### S CYCLES






- **134°C Universal Cycle:** 3,5 min. of sterilisation plus the drying phase
- **121°C Softprogramm Cycle:** 15,5 min. of sterilisation plus the drying phase

#### SPECIAL CYCLES



- **Drying cycles:** drying phase only

NOTE: For each available sterilisation cycles, the manufacturer must perform the Works Tests in accordance with EN 13060 (par. 7.3). Any available cycles for which the manufacturer has not performed the Works Tests must undergo Installation Tests in accordance with EN 13060 par. 7.4.

### 4.3.2 Diagrammi Cicli di Tipo B

### 4.3.3 Cycle Diagram type S

## 4.4 Sterilisation program specifications

|                                      |            | Sterilisation cycles type B |            |         |        |          | Sterilisation cycles type S |              |
|--------------------------------------|------------|-----------------------------|------------|---------|--------|----------|-----------------------------|--------------|
| Cycle name                           |            | Standard                    | Standard + | B-Flash | Prion  | Standard | Universal S                 | Softprogramm |
| Sterilisation temperature            |            | 134°C                       |            |         |        | 121°C    | 134°C                       | 121°C        |
| Sterilisation pressure               |            | 316 kPa                     |            |         |        | 214kPa   | 316 kPa                     | 214kPa       |
| Number of vacuum phases              |            | 3                           | 4          | 3       | 4      | 3        | 2                           | 2            |
| Sterilisation phase duration         |            | 5,5 min                     | 5,5 min    | 5,5 min | 20 min | 20,5 min | 3,5 min                     | 15,5 min     |
| Minimum drying phase duration        |            | 16 min                      | 16 min     | 6 min   | 16 min | 16 min   | 16 min                      | 16 min       |
| Average total cycle duration         | HC 50      | 45 min                      | 49 min     | 33 min  | 60 min | 58 min   | 36 min                      | 47 min       |
|                                      | HC 50 Plus | 48 min                      | 51 min     | 34 min  | 62 min | 60 min   | 37 min                      | 48 min       |
| Average H <sub>2</sub> O consumption | HC 50      | 415 ml                      | 430 ml     | 300 ml  | 525 ml | 455 ml   | 260 ml                      | 295 ml       |
|                                      | HC 50 Plus | 455 ml                      | 470 ml     | 315 ml  | 565 ml | 500 ml   | 275 ml                      | 300 ml       |
| Average energy consumption           | HC 50      | 515 Wh                      | 530 Wh     | 390 Wh  | 600 Wh | 547 Wh   | 375 Wh                      | 410 Wh       |
|                                      | HC 50 Plus | 646 Wh                      | 655 Wh     | 425 Wh  | 710 Wh | 654 Wh   | 450 Wh                      | 480 Wh       |

## 4.5 Putting the material to be sterilised in the chamber

Package the load very carefully, making sure that the medical devices are compatible with the chosen packaging.

The chosen packaging material must be:

- Compatible with the sterilisation process.
- Compatible with the labelling system.

Only use medical packaging that is compliant with: EN 868-5 and ISO 11607-1 qualified for 134°C cycles.

To check that the thermodynamic sterilisation cycle has completed successfully, always put a class 5 chemical indicator in together with the load in accordance with ISO 11140-1. Choose the indicator according to the selected cycle (121°C or 134°C)

Place the clean material on the trays and distribute them evenly.

The chamber can contain a maximum of:

|            | B Cycles                             |                       |             |                       |             | S Cycles              |             |
|------------|--------------------------------------|-----------------------|-------------|-----------------------|-------------|-----------------------|-------------|
|            | Standard - Prion - Standard + Cycles |                       |             | 134 B-Flash Cycles    |             |                       |             |
|            | Unpacked solid material              | Packed solid material | Porous load | Packed solid material | Porous load | Packed solid material | Porous load |
| HC 50      | 5,5 kg                               | 3,5 kg                | 1,5 kg      | 0,5 kg                | 0,2 kg      | 2,5 kg                | 0,3 kg      |
| HC 50 Plus | 6,5 kg                               | 4 kg                  | 2 kg        | 0,7 kg                | 0,3 kg      | 3,0 kg                | 0,4 kg      |

Objects in bags must be placed on the trays with the transparent part of the package facing down.

### 4.5.1 Preparing the material

First of all, remember that when handling and moving contaminated material, it is good practice to take the following precautions:

- Wear suitably thick rubber gloves and a face mask;
- Wash your already-gloved hands with a germicidal detergent;
- Always use a tray to transport instruments.
- Never transport then directly in your hand;
- To avoid the risk of contracting hazardous infections, protect you hands from contact with any sharp or pointed parts;
- Immediately separate any items that must not be sterilised or that cannot support the process;
- When you have finished handling the unsterile material, wash your still-gloved hands thoroughly. All materials and/or instruments to be sterilised must be perfectly clean and free of any kind of residue (organic/inorganic deposits, paper fragments, cotton/gauze pads, limescale, etc.).

#### NOTE

In addition to causing problems during the sterilisation process, failure to clean and remove residues can damage the instruments and/or the steriliser.

#### NOTE

In addition to causing problems during the sterilisation process, failure to clean and remove residues can damage the instruments and/or the steriliser.

Proceed as described to clean effectively:

1. Rinse the instruments under a jet of running water immediately after use;
2. Separate metal instruments by material type (carbon steel, stainless steel, brass, aluminium, chromium, etc.) to prevent electrolytic oxidation;
3. Wash with an ultrasonic cleaner containing a mixture of water and germicidal solution, carefully following the manufacturers recommendations. Alternatively, use a thermal disinfectant.
4. For better results, use a detergent specifically for ultrasonic cleaning, with a neutral pH.

#### NOTE

Solutions that contain phenol or quaternary ammonium based compounds can corrode the instruments and the metal parts of the ultrasonic equipment.

5. After washing, rinse the instruments thoroughly and check that the residues have been completely removed. Repeat the washing cycle or wash by hand if necessary.

#### NOTE

To prevent the formation of limescale stains, rinse with deionised or distilled water if possible. Always dry the instruments if you use hard tap water.

For handpieces (turbines, counterangles, etc.), in addition to the above, use special equipment that will clean them properly inside (which may also include lubrication).

#### NOTE

When the sterilisation program has finished, remember to lubricate the handpiece internal mechanisms with special oil. These precautions will ensure that the lifetime of the instrument is not reduced.

#### CAUTION

Refer to the instrument/material manufacturer's instructions and make sure that it is compatible before treating it in the autoclave. Strictly follow the usage methods for detergents and disinfectants and the operating instructions for automatic washing and/or lubrication equipment.

For textiles, (or porous materials in general), such as lab coats, napkins, caps and similar, wash them thoroughly and dry them before treating them in the autoclave.

#### NOTE

Do not use strong detergents that contain chlorine and/or phosphates. Do not bleach with chlorine-based products. These components could damage the tray supports, the trays and any metal instruments in the sterilisation chamber.

### 4.6 Filling and draining water

Volume of water used for a sterilisation cycle/load with the highest consumption: 700ml

Minimum water tank filling level: 1.1 litres

#### 4.6.1 Filling the feed tank

- Only use water that complies with the technical specifications given in 7.1 and the requirements according to EN 13060 Appendix C.
- The display shows MSG001 (empty load tank), fill the feed tank in one of the following two ways:

##### Automatic filling from the front

- Insert the part of the hose with the quick coupling into the water filling fitting (6 fig. 1) and put the other end into the water container.
- Select the "water filling" menu .
- Press the icon "PLAY WATER FILLING"  to start the water filling pump.
- Filling stops automatically when the water reaches the maximum level in the feed tank .
- Press the icon "STOP"  to voluntarily stop filling.

##### REAR AUTOMATIC FILLING

In case of connection to the demineralization system (18 fig. 2) the filling of the charging tank is automatic at the beginning and end of the sterilization cycle. During the water loading the machine cannot perform sterilization cycles and tests.

- 01 Select the "Water filling" menu
- 02 Select "Start" to start filling
- 03 Select "Stop" to stop filling
- 04 The cursor position indicates the water quality

N.B. Filling will stop automatically when the maximum level has been reached, and you will return to the Home screen

#### Manual loading tank filling

Open the cover of the loading tank identified.  
Fill the tank with water according to the specifications provided. Do not exceed the maximum level line on the tank.  
Close the cover.

#### Emergency manual filling

- Insert the part of the hose with the quick coupling into the hole (9 fig. 1)
- Insert the funnel into the other end of the hose and pour in a maximum of 2 litres of distilled water.
- In case of connection to demineralization system (18 fig. 2) the loading tank is filled automatically.

#### 4.6.2 Draining the discharge tank

- Insert the end of the hose without the quick coupling into a container to collect the drained water.
- Insert the quick coupling into the hole (7 fig. 1) and position the container under the autoclave to drain the tank.
- Drain the water in compliance with the waste disposal regulations in force.
- Possibility to use the drain tank outlet, placed at the back (15 fig. 2).

## 4.7 Test programs

- 01 Select the "Test cycles" menu
- 02 Alternatively, directly select "Last test run" and go to step 04
- 03 Select the desired cycle
- 04 Screen with all the parameters of the selected cycle
- 05 (Optional) you can open the door to insert other instruments
- 06 (Optional) Select the programmed start button to set the start time
- 07 Start the test cycle
- 08 (If anabled) enter the user PIN.
- On confirming, the test cycle will start.

- 09 You can use the green bar to monitor the cycle progress.
- 10 When the cycle finishes, press the confirm icon to return to the Home screen.
- 11-12 You can stop the cycle manually by holding down the "STOP CYCLE" button for 3 seconds. Wait for the machine to be put in a safe condition and press the confirm icon (12) and then press Home to return to the home screen.
- 14/15 Cycle Validation Protocol: The selection of the Batch, Load and/or Indicator validation options generates the relative information in the cycle report: "Batch Ok" – "Batch Fail" "Indicator OK" – "Indicator Fail" – "No Indicator". "Load OK" – "Load Fail" – "No Load" \* The confirmation button (15) will be "Operational" when the user makes a selection.
- 16 For the validation of Batch, Load and/or Indicator, the user must enter his PIN. When the confirmation button is pressed (16), the system generates the Cycle Report PDF with the identification of the User who performed the cycle validation.

### 4.7.1 Test program types

Selecting the “test” sub-menu accesses the list of test cycles that autoclave can run:

- Vacuum Test
- 134°C Helix Test
- 121°C Helix Test
- 134°C Bowie&Dick Test

### 4.7.2 Vacuum test (leakage test)

To run this function, the chamber must be cold and have no instruments in it.




### 4.7.3 Bowie & dick test

To run this test, the chamber must have no instruments in it. The Bowie & Dick “package” used for the porous load test must be positioned on the lower tray.

### 4.7.4 Helix test

To run this test, the chamber must have no instruments in it. The instrument used to test narrow hollow loads must be positioned on the lower tray.

### 4.7.5 Test diagrams


## 4.8 Settings


### 4.8.1 User settings

- 01 Select "Settings"
- 02 Select "User"
- 03 Enter the administrator pin (factory setting "2222") and confirm  
With the Administrator PIN (2222), it is possible to set the users and to manage their PIN and the Cycle Validation Protocol.
- 04 Click on the arrows to enable or disable the USERS function.  
When the function is enabled, a user PIN will be requested every time a cycle/ test is started.  
If PIN User Management is enabled, upon confirmation you will go to the new Cycle Validation Protocol Management screen. If PIN User Management is disabled, upon confirmation of the selection, the confirmation button ④ will become "Non-Operational".
- 05 When the "Users" function is enabled, you can ADD other users
- 06 When the "Users" function is enabled, you can DELETE previously created users
- 07 When the "Users" function is enabled, you can EDIT previously created users
- 08 Confirm to save
- 09 By enabling the Protocol Management PIN, the Cycle Validation Protocol is activated at the end of the cycle.  
The Confirmation button ⑨ will be "Operational" when the user changes the Disable/Enable state.

### 4.8.2 Date and time settings

- 01 Select "Settings"
- 02 Select "date and time"
- 03 Click on the arrows to increase or decrease the relative field (day, month, year)
- 04 Confirm to go to the "time" screen
- 05 Click on the arrows to increase or decrease the relative field (hours, minutes, 12/24)
- 06 Confirm to save


#### 4.8.3 Printer settings

- 01 Select "Settings"
- 02 Select "Printer"
- 03a Click on the arrows to change from automatic to manual print mode and vice versa
- 03b Click on the arrows to select the type of the printing (Report – Label – Barcode Label)
- 03c Click on the arrows to select the Model of the printer (Hygoprint Easy – MCP7810X)
- 04 Confirm to save


#### 4.8.4 Language settings

EXAMPLE OF  
LANGUAGES  
SCREEN

- 01 Select "Settings"
- 02 Select "Language"
- 03 Click on the desired language

#### 4.8.5 Measurement unit settings


- 01 Select "Settings"
- 02 Select "Measurement unit"
- 03 Click on the arrows to change from "bar" to "kPa" to "mbar"
- 04 Click on the arrows to change from "absolute" to "relative"
- 05 Confirm to save

#### 4.8.6 Water settings

- 01 Select "Settings"
- 02 Select "Water settings"
- 03 Click on the arrows to change from front filling (manual) to rear filling (from the demineralization system)
- 04 Confirm to save

## ETHERNET

- 01 Select "Settings"
- 02 Select "Connectivity"
- 07 Select "ETHERNET"
- 08/09/10 Click on the relative fields and enter the data
- 11 Confirm to save

## 4.8.8 Hygoclave 50 / Hygoclave 50 Plus - Ethernet Connection

The Hygoclave 50/Hygoclave 50 Plus series sterilizers can be connected to a local Ethernet network by means of a specific cable connected to the special connector (14 - fig. 2) located in the rear of the autoclave.


Through a PC or other device connected to the network (smartphone, tablet, etc.) with the use of a Web Browser (Internet browsing program like Internet Explorer, Firefox, Chrome, etc.) and by knowing the TCP-IP number assigned to the sterilizer, you will be able to monitor the machine to know its operating status and download the reports of the cycles performed.

## ETHERNET NETWORK CONNECTION CONFIGURATION

The sterilizer Web Server only works with a network with active DHCP and static IP.

The number chosen as the machine address must be compatible with the range assigned by the DHCP server of the local network and its Subnet Mask.

To find out the local network numbering range on Windows systems, you can use the IPCONFIG command from a "Command Prompt" window (accessible from "All Programs – Accessories"):


For a correct configuration you must assign to the IP address of the machine: - a number belonging to the local network defined by the first 3 values of the IP address: in the example above it is "169.254.190.xxx"; - as the last number a value not assigned to the local network: e.g. "169.254.190.158". Therefore the complete IP address of the sterilizer could be of the type: **169.254.190.158**

Once you know the IP address and Subnet Mask to be assigned to the machine you must enter these values in the setting screen of the autoclave Ethernet parameters using the keyboard and save the data:


Note: The Gateway address is not important for local network communications.

Once the Ethernet configuration data has been entered and saved, the autoclave must be switched off and on again so that the data can be used by the Web Server of the machine.



At this point, with the sterilizer turned on on the Home screen, from a PC connected to the local network, or directly connected to the machine via Ethernet cable, you must open the Web Browser (Internet browsing program like Internet Explorer, Firefox, Chrome, etc.) and in the address field you must enter the IP address previously entered in the sterilizer (in our e.g. **169.254.190.158**). You will then have access to a "Hygoclave 50 / Hygoclave 50 Plus control panel" web page dedicated to downloading reports of cycles that the machine has performed (**Download**) or to monitor its operation (**Monitor**):

Hygoclave 50 / Hygoclave 50 Plus control panel


Hygoclave 50 / Hygoclave 50 Plus report download page

Hygoclave 50 / Hygoclave 50 Plus monitor page

## 4.9 Data output

### USB


- 01 Select the “Data output” menu
- 02 Select the USB symbol
- 03 Insert the USB pen drive and wait for it to be recognised
- 04 Select the number of cycles to transfer
- 05 Confirm and wait for the transfer to complete
- 13 Download progress bar-graph

### PRINTER

- 01 Select the “Data output” menu
- 06 Select the PRINTER symbol
- 07/08/09 Select the printout type (report, labels, barcodes)
- 10 Select the number of labels/barcodes to print
- 11 Confirm printout
- 12 Setting of the Expiry Date

### INFO

- 01 Select the “Data output” menu
- 05 Select the INFO symbol

## DISPLAY CALIBRATION

- [illegible]

EXAMPLE PRINTOUT

EN

|                                                                                                                    |  |                       |                     |                                                                                                                                                           |
|--------------------------------------------------------------------------------------------------------------------|--|-----------------------|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>DÜRR DENTAL SE</b><br>Höpfheimer Str. 17, 74321<br>Bietigheim-Bissingen<br>Phone +49 07142 705 0                |  |                       |                     |                                                                                                                                                           |
| <b>Firmware</b><br>firmware version                                                                                |  |                       |                     |                                                                                                                                                           |
| <b>Serial Number</b><br>aa m STN xxx nnnnnn                                                                        |  |                       |                     | aa= year; m= month; xxx= model;<br>nnnnn= progressive no.                                                                                                 |
| <b>Water Quality</b><br>conductibility                                                                             |  |                       |                     |                                                                                                                                                           |
| <b>Program</b><br>134 Standard                                                                                     |  |                       |                     | program name                                                                                                                                              |
| <b>Start</b><br>Date      Time                                                                                     |  |                       |                     |                                                                                                                                                           |
| <b>Time</b><br>(hh:mm:ss)                                                                                          |  | <b>Press</b><br>(kPa) | <b>Temp</b><br>(°C) |                                                                                                                                                           |
| <b>Pre-heating</b><br>time      pressure      temperature                                                          |  |                       |                     | values at the end of the preheating phase                                                                                                                 |
| <b>Fractional Vacuum</b><br>time      pressure      temperature                                                    |  |                       |                     | values on reaching the 1st vacuum threshold                                                                                                               |
| time      pressure      temperature                                                                                |  |                       |                     | values on reaching the pressure threshold                                                                                                                 |
| np = xxx                                                                                                           |  |                       |                     | np= xxx number of injection pump pulses                                                                                                                   |
| time      pressure      temperature                                                                                |  |                       |                     | values on reaching the 2nd vacuum threshold                                                                                                               |
| time      pressure      temperature                                                                                |  |                       |                     | values on reaching the pressure threshold                                                                                                                 |
| np = xxx                                                                                                           |  |                       |                     | np= xxx number of injection pump pulses                                                                                                                   |
| time      pressure      temperature                                                                                |  |                       |                     | values on reaching the 3rd vacuum threshold                                                                                                               |
| time      pressure      temperature                                                                                |  |                       |                     | values on reaching the pressure threshold                                                                                                                 |
| np = xxx                                                                                                           |  |                       |                     | np= xxx number of injection pump pulses                                                                                                                   |
| <b>Equilibration</b><br>time      pressure      temperature                                                        |  |                       |                     | initial values                                                                                                                                            |
| time      pressure      temperature                                                                                |  |                       |                     | final values                                                                                                                                              |
| <b>Sterilization</b><br>time      pressure      temperature                                                        |  |                       |                     | print values every xx seconds (with<br>xx equal to 10s or 15s or 30s or 60s)                                                                              |
| ...      ...      ...                                                                                              |  |                       |                     |                                                                                                                                                           |
| ...      ...      ...                                                                                              |  |                       |                     |                                                                                                                                                           |
| ...      ...      ...                                                                                              |  |                       |                     |                                                                                                                                                           |
| ...      ...      ...                                                                                              |  |                       |                     |                                                                                                                                                           |
| time      pressure      temperature                                                                                |  |                       |                     | final values                                                                                                                                              |
| <b>Pmax</b> =      ...                                                                                             |  |                       |                     | maximum and minimum values                                                                                                                                |
| <b>Pmin</b> =      ...                                                                                             |  |                       |                     | during the entire sterilisation                                                                                                                           |
| <b>Tmax</b> =      ...                                                                                             |  |                       |                     |                                                                                                                                                           |
| <b>Tmin</b> =      ...                                                                                             |  |                       |                     |                                                                                                                                                           |
| <b>np</b> =      ...                                                                                               |  |                       |                     | number of injection pump pulses<br>during the sterilisation                                                                                               |
| <b>Discharge</b><br>time      pressure      temperature                                                            |  |                       |                     | value at the end of the phase                                                                                                                             |
| <b>Dry</b><br>time      pressure      temperature                                                                  |  |                       |                     | value at the end of the phase                                                                                                                             |
| npv = yyy                                                                                                          |  |                       |                     | npv= yyy number of vacuum pulses                                                                                                                          |
| <b>Pressure Balance</b><br>time      pressure      temperature                                                     |  |                       |                     | value at the end of the phase                                                                                                                             |
| <b>End</b><br>Date      Time                                                                                       |  |                       |                     |                                                                                                                                                           |
| Duration                                                                                                           |  |                       |                     | duration of the entire cycle                                                                                                                              |
| <b>OK</b>                                                                                                          |  |                       |                     |                                                                                                                                                           |
| <b>Cycle Number</b><br>aaaaa/bbbbb (ccccc)                                                                         |  |                       |                     | aaaaa= no. of successful cycles;<br>bbbbbb= no. of cycles started (total);<br>ccccc= no. of cycles started for the<br>specific cycle                      |
| <b>User</b><br>Administrator      (Name and Surname)                                                               |  |                       |                     | user name of user management if active                                                                                                                    |
| <b>Signature</b><br><br>-----                                                                                      |  |                       |                     |                                                                                                                                                           |
| <b>Validation protocol</b><br>Batch OK or Fail<br>Load OK or Fail - No Load<br>Indicator OK or Fail - No Indicator |  |                       |                     | Information relating to Batch, Load<br>and Indicator validation. This data is<br>present only if the user has activated<br>the Cycle Validation Protocol. |
| <b>User</b><br>Administrator      (Name and Surname)                                                               |  |                       |                     | Name of the user that validated the<br>cycle data (if the Cycle Validation<br>Protocol has been activated).                                               |
| <b>Signature</b><br><br>-----                                                                                      |  |                       |                     |                                                                                                                                                           |

## 5 ROUTINE MAINTENANCE

### 5.1 General table

| Operation                                               | Frequency         |                      |                         |                         |
|---------------------------------------------------------|-------------------|----------------------|-------------------------|-------------------------|
|                                                         | Daily<br>5 cycles | Monthly<br>50 cycles | Quarterly<br>500 cycles | Annually<br>1000 cycles |
| Clean the door seal                                     | •                 |                      |                         |                         |
| Clean the chamber edge<br>(area where the seal touches) | •                 |                      |                         |                         |
| Clean the door cover                                    |                   | •                    |                         |                         |
| Clean the chamber and accessories                       |                   | •                    |                         |                         |
| Clean the front panel                                   |                   | •                    |                         |                         |
| Clean the side walls and top panel                      |                   | •                    |                         |                         |
| Clean the chamber filter                                |                   | •                    |                         |                         |
| Clean the dust filter (if provided)                     |                   | •                    |                         |                         |
| Clean the screw and nut screw                           |                   | •                    |                         |                         |
| Lubricate the screw and nut screw                       |                   | •                    |                         |                         |
| Clean the feed tank                                     |                   |                      | •                       |                         |
| Replace the bacterial filter                            |                   |                      | •                       |                         |
| Replace the door seal                                   |                   |                      |                         | •                       |

Before cleaning any surface, make sure that it has cooled down.

#### 5.1.1 Cleaning the sterilisation chamber and internal components

Use a soft cloth soaked in water to avoid damaging the protective surface treatment of the chamber. Rinse thoroughly with water and then dry. If necessary, moisten the cloth with white vinegar to remove any limescale.

However, we recommend cleaning the chamber only when strictly necessary. If there are limescale deposits, immediately check what is causing them. White patches forming at the bottom of the chamber means that the water used is of poor quality or does not comply with the required specifications.

#### 5.1.2 Cleaning the door cover

Use a soft cloth soaked in water to avoid damaging the protective treatment of the door cover. Rinse thoroughly with water and then dry. **DO NOT USE** metallic or particularly abrasive cloths for cleaning. Using them could damage the external surface of the door cover.

#### 5.1.3 Cleaning the external front panel

Use a soft cloth soaked in alcohol, taking care not to let it penetrate into the autoclave internal parts. Never use detergents based on triline, benzene or turpentine or solvents of any kind. **DO NOT USE** metallic or particularly abrasive cloths for cleaning.


#### 5.1.4 Cleaning the door seal

Clean the seal with a soft cloth moistened with water or white vinegar, to remove any limescale. Rinse thoroughly if you used vinegar. This procedure serves to remove impurities that can cause pressure loss in the sterilisation chamber and possibly break the seal.

#### 5.1.5 Cleaning the filter in the chamber

- Remove the tray holder with the trays from the chamber
- Rinse the filter thoroughly with water.
- Reposition the filter on the bottom of the chamber and screw it back in place.

If the cleaning operation does not produce satisfactory results, call DÜRR DENTAL Technical Support.

#### 5.1.6 Replacing the Bacterial or Bio X filter

Open the autoclave door

- Unscrew the filter (**2 Fig. 1**) by turning it anticlockwise.
- Screw the new filter in clockwise, making sure that it is tightened properly.

#### 5.1.7 Replacing the door seal



This must be carried out by qualified technicians.

#### 5.1.8 Cleaning the dust filter (if provided)

Remove the dust filter from the bottom of the autoclave (**19 fig. 1**), rinse it thoroughly with water and dry it before refitting it.

### 5.1.9 Lubricating and cleaning the screw and nut screw

Every 50 cycles, clean and lubricate with silicone grease the inside of the threaded bushing "A" located on the door .

Clean screw "B" on the front panel using a soft cloth.

⊘ Never use hydrochloric acid or detergents containing chlorine to clean any part of the autoclave.

### 5.1.10 Drain tank cleaning and disinfection

- Carry out the following steps:
- Empty the drain and loading tanks.
- Remove all the connections fastened to the loading and drain tanks.
- Disinfect the inside surfaces of the water tanks by wiping them with a surface disinfectant, e.g. Dürr Dental FD 322 or FD 333.
- Once the inside surfaces of the tanks have dried, fill the loading tank with demineralized water (see section 7.1).
- Run 2 fast program cycles.
- Empty the loading tank and dry it with a sterile cloth.

### 5.2 Maintenance

After every 1000 cycles or two years, a reminder message appears to perform maintenance.

The maintenance must only be carried out by an authorised and qualified service technician. The maintenance consists of replacing, cleaning and checking some parts, as indicated in the relative table.

The work consists of the following:

#### **Replace:**

- door seal
- bacterial or Bio X filter
- chamber filter
- water filter

#### **Check:**

- Safety valve check
- Leakage current check
- Hose checks: hoses maintain performance until they physically break, in which case they do not work. Only replace them with DÜRR DENTAL spare parts.

When finished, the qualified technician will issue a certificate stating that the autoclave maintenance was carried out.

Comply with the intervals specified or indicated in this manual and in the additional sheets. Only use ORIGINAL SPARE PARTS. Failure to comply will void the equipment warranty.

### 5.3 Periodic process effectiveness and safety checks

| Operation                      | Frequency        |       |          |                 |
|--------------------------------|------------------|-------|----------|-----------------|
|                                | Before first use | Daily | Annually | Every two years |
| Vacuum Test EN 17665           |                  | •     |          |                 |
| Bowie & Dick Test 17665        |                  | •     |          |                 |
| Helix Test EN 17665            |                  | •     |          |                 |
| Protective earthing continuity |                  |       |          | •               |
| Insulation test                |                  |       |          | •               |
| Casing leakage current         |                  |       |          | •               |
| Initial validation EN 17665    | •                |       |          |                 |
| Periodic validationa EN 17665  |                  |       | •        |                 |
| Safety valve check             |                  |       | •        |                 |

Before proceeding with the operations indicated above, always switch off the main switch (5 Fig. 1) and make sure that the chamber is cold.

Never use hydrochloric acid or detergents containing chlorine to clean any part of the autoclave.

## 6 Troubleshooting guide

### 6.1 Error Codes

| Code   | Message                                    | Description                                                                                  | Reset procedure |
|--------|--------------------------------------------|----------------------------------------------------------------------------------------------|-----------------|
| ER 001 | Door closing timeout                       | The door closing time is too long                                                            | B               |
| ER 002 | Door switch                                | Door switch error condition                                                                  | B               |
| ER 003 | Door unlocking timeout                     | The door unlocking time is too long                                                          | B               |
| ER 004 | Door opening timeout                       | The door opening time is too long                                                            | B               |
| ER 006 | Very poor water quality                    | The water conductivity is higher than 30uS/cm                                                | B               |
| ER 007 | Floating feed tank in unexpected condition | Floating feed tank in unexpected condition                                                   | B               |
| ER 100 | Vacuum timeout during preheating           | The vacuum preheating phase time is too long                                                 | A               |
| ER 101 | Pressure rise timeout during preheating    | The vacuum preheating phase time is too long                                                 | A               |
| ER 110 | Vacuum pressure timeout                    | The time to reach the vacuum pressure is too long                                            | A               |
| ER 120 | Discharge pressure timeout                 | The time to reach the thermodynamic discharge pressure is too long                           | A               |
| ER 121 | Pressure difference                        | The pressure increase is lower than the set parameter                                        | A               |
| ER 130 | Discharge time                             | The pressure discharge time is too long                                                      | A               |
| ER 140 | Sterilisation pressure timeout             | The time to reach the sterilisation pressure is too long                                     | A               |
| ER 141 | Pressure difference                        | The pressure increase is lower than the set parameter                                        | A               |
| ER 150 | T chamber high                             | The chamber temperature is over the top limit for sterilisation                              | A               |
| ER 151 | T chamber low                              | The chamber temperature is below the bottom limit for sterilisation                          | A               |
| ER 152 | T theoretical high                         | The theoretical temperature is over the top limit for sterilisation                          | A               |
| ER 153 | T theoretical low                          | The theoretical temperature is below the bottom limit for sterilisation                      | A               |
| ER 154 | Temperature difference                     | The difference between the chamber and theoretical temperature is over the allowed limit     | A               |
| ER 160 | Final discharge                            | The final discharge time is too long                                                         | A               |
| ER 171 | Vacuum pressure timeout                    | The time to reach the vacuum pressure is too long                                            | A               |
| ER 172 | AV block                                   | The measured pressure is below the set parameter for the drying phase                        | A               |
| ER 173 | Vacuum pressure timeout                    | The time to reach the vacuum pressure is too long                                            | A               |
| ER 180 | Balancing timeout                          | The pressure balancing time is too long                                                      | A               |
| ER 200 | Vacuum test timeout                        | The time to reach the vacuum pressure is too long in the vacuum test                         | A               |
| ER 201 | Vacuum test: first phase                   | Failure to comply with the parameters in the first phase of the vacuum test                  | A               |
| ER 202 | Vacuum test: second phase                  | Failure to comply with the parameters in the 2nd phase of the vacuum test                    | A               |
| ER 203 | Vacuum test: temperature                   | Failure to comply with the temperature limits during the vacuum test                         | A               |
| ER 204 | Vacuum test: balancing                     | Pressure balancing failed in the vacuum test                                                 | A               |
| ER 310 | Vacuum pressure threshold reached timeout  | The time to reach the test vacuum pressure threshold has been exceeded                       | A               |
| ER 340 | Test pressure rise timeout                 | The time to reach the test vacuum pressure threshold has been exceeded                       | A               |
| ER 341 | Pressure delta                             | Pressure delta                                                                               | A               |
| ER 342 | Safety valve open                          | The safety valve has opened during test execution                                            | A               |
| ER 355 | Test execution timeout                     | The test execution time has been exceeded                                                    | A               |
| ER 360 | Final discharge timeout                    | The time to discharge the pressure has been exceeded                                         | A               |
| ER 370 | AV lock during drying                      | Pressure measurement less than the parameter established in the drying phase during the test | A               |
| ER 380 | Pressure balancing timeout                 | The pressure balancing time has been exceeded during test execution                          | A               |

## 6.2 Alarm Codes

| Code   | Message                                              | Description                                                                        | Reset procedure |
|--------|------------------------------------------------------|------------------------------------------------------------------------------------|-----------------|
| AL 001 | Reset pressed                                        | The reset button was pressed                                                       | A/B             |
| AL 002 | Blackout                                             | A blackout has occurred                                                            | A/B             |
| AL 003 | Water filling time                                   | The time to fill the feed tank is too long                                         | A               |
| AL 100 | Chamber probe fault                                  | The chamber probe reading circuit is faulty                                        | A               |
| AL 101 | Generator probe fault                                | The generator probe reading circuit is faulty                                      | A               |
| AL 102 | Bundle probe fault                                   | The pipe bundle probe reading circuit is faulty                                    | A               |
| AL 103 | Pressure transducer                                  | The pressure transducer is faulty                                                  | A               |
| AL 104 | Door microswitch                                     | The door microswitch has opened unexpectedly                                       | A               |
| AL 105 | Safety microswitch                                   | The safety microswitch has opened unexpectedly                                     | A               |
| AL 200 | Chamber safety cut-off                               | The chamber has exceeded the safety temperature                                    | A               |
| AL 201 | Generator safety cut-off                             | The generator has exceeded the safety temperature                                  | A               |
| AL 202 | Bundle safety cut-off                                | The pipe bundle has exceeded the safety temperature                                | A               |
| AL 203 | Safety valve fault                                   | The pressure has exceeded the maximum trip threshold of the safety valve           | A               |
| AL 300 | Hardware error<br>Reboot the system                  | There is no synchronisation for managing AC loads and checking the mains frequency | A               |
| AL 301 | Valve drive voltage                                  | There is no valve drive voltage due to a PTC trip                                  | A               |
| AL 302 | Mains frequency error                                | The mains frequency is < 45 Hz or > 65 Hz                                          | A               |
| AL 303 | Watchdog trip                                        | Program execution has crashed                                                      | A               |
| AL 304 | Communication error with Main Board<br>Reboot system | There is a communication fault between the main board and the CPU board            | A               |
| AL 305 | Hardware error                                       | Hardware fault on the main board                                                   | A               |
| AL 306 | Hardware error                                       | Hardware fault on the main board                                                   | A               |
| AL 307 | Hardware error                                       | Hardware fault on the main board                                                   | A               |
| AL 308 | Hardware error                                       | Hardware fault on the main board                                                   | A               |
| AL 309 | Hardware error                                       | Hardware fault on the main board                                                   | A               |
| AL 310 | Hardware error                                       | Hardware fault on the main board                                                   | A               |
| AL 311 | Hardware error                                       | The cycle has incorrectly been closed – Software error                             | A               |
| AL 401 | Generator temperature low                            | Generator failure                                                                  | A               |
| AL 402 | Band temperature low                                 | Band failure                                                                       | A               |

## 6.3 Messages

| Code   | Phase                                   | Message                                                   |
|--------|-----------------------------------------|-----------------------------------------------------------|
| Msg001 | Standby or cycle start                  | Feed tank empty                                           |
| Msg002 | Standby or cycle start                  | Feed tank full                                            |
| Msg003 | Pin entry                               | Incorrect Pin                                             |
| Msg004 | Water loading                           | Poor water quality                                        |
| Msg005 | Cycle start                             | Door open                                                 |
| Msg008 | Vacuum Test start                       | Chamber temperature over limit                            |
| Msg009 | Machine switch-on                       | Serial number inconsistency between CPU and MB            |
| Msg010 | Machine switch-on                       | Lubricate the automated system of the door                |
| Msg011 | Machine switch-on                       | Replace the Bio-X filter                                  |
| Msg012 | Machine switch-on                       | Conduct maintenance                                       |
| Msg013 | Incorrect date                          | System date incorrect - Check                             |
| Msg014 | Do vacuum test                          | Do vacuum test                                            |
| Msg016 | PIN too short                           | PIN too short                                             |
| Msg017 | Duplicate PIN                           | PIN already exists                                        |
| Msg018 | MB FW revision incompatible with CPU FW | Firmware version incompatible – contact Technical Support |

## 6.4 Reset procedures

### Reset procedure "A"

If the cycle is "running", reset the error by selecting the confirm icon on the safety screen and then select the Home icon to return to the initial screen.

### Reset procedure "B"

Reset the error by selecting the confirm icon on the safety screen.

**N.B. In the event of faults that compromise sterilisation/test cycle execution, the cycle selection icons on the Home screen will be disabled.**

## 6.5 Troubleshooting

| Code   | Possible cause                                                                                   | Possible solution                                                                                                                                                                                                                                                      |
|--------|--------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ER 001 | There is a malfunction in the motorised door closing/opening system.                             | Reset the machine.<br>Repeat the door closing procedure.<br>If the problem persists, contact technical support.                                                                                                                                                        |
| ER 002 | There is a malfunction in the motorised door closing/opening system.                             |                                                                                                                                                                                                                                                                        |
| ER 003 | There is a malfunction in the motorised door closing/opening system.                             |                                                                                                                                                                                                                                                                        |
| ER 004 | There is a malfunction in the motorised door closing/opening system.                             |                                                                                                                                                                                                                                                                        |
| ER 006 | The water quality is very poor.                                                                  | Reset the machine.<br>Very poor water quality: empty the feed tank and refill it with better quality water.<br>If the problem persists, contact technical support.                                                                                                     |
| ER 007 | Floating feed tank in unexpected condition.                                                      | Reset the machine. Check the condition of the floating feed tank.<br>If the problem persists, contact Technical Support.                                                                                                                                               |
| ER 100 | There is a leak in the machine plumbing circuit or a malfunction in the vacuum pump.             | Reset the machine. Check that the door seal is clean and positioned correctly, then repeat the selected program. If the problem persists, contact technical support.                                                                                                   |
| ER 101 | There is a leak in the machine plumbing circuit.                                                 | Reset the machine. Check:<br>a) that the door seal is clean and positioned correctly;<br>b) that there is water in the feed tank;<br>c) the amount of load in the chamber.<br>Then repeat the selected program.<br>If the problem persists, contact technical support. |
| ER 110 | There is a leak in the machine plumbing circuit or a malfunction in the vacuum pump.             | Reset the machine. Check:<br>a) that the door seal is clean and positioned correctly;<br>b) that the filter in the chamber is clean and positioned correctly.<br>Then repeat the selected program.<br>If the problem persists, contact technical support.              |
| ER 120 | There is a leak in the machine plumbing circuit or a malfunction in the water injection pump.    | Reset the machine. Check:<br>a) that the door seal is clean and positioned correctly;<br>b) that there is water in the feed tank;<br>c) the amount of load in the chamber.<br>Then repeat the selected program.<br>If the problem persists, contact technical support. |
| ER 121 | There is a leak in the machine plumbing circuit or a malfunction in the water injection pump.    |                                                                                                                                                                                                                                                                        |
| ER 130 | There is a malfunction in the machine discharge system.                                          | Reset the machine.<br>Check the amount of load in the chamber and repeat the selected program.<br>If the problem persists, contact technical support.                                                                                                                  |
| ER 140 | There is a leak in the machine plumbing circuit or a malfunction in the water injection pump.    | Reset the machine. Check:<br>a) that the door seal is clean and positioned correctly;<br>b) that there is water in the feed tank;<br>c) the amount of load in the chamber.<br>Then repeat the selected program.<br>If the problem persists, contact technical support. |
| ER 141 | There is a leak in the machine plumbing circuit or a malfunction in the water injection pump.    |                                                                                                                                                                                                                                                                        |
| ER 150 | There is a leak in the machine discharge system or a chamber loading error.                      | Reset the machine. Check:<br>a) the tray support position and the tray position in the chamber;<br>c) the amount of load in the chamber.<br>Then repeat the selected program.<br>If the problem persists, contact technical support.                                   |
| ER 151 | There is a malfunction in the sterilisation condition holding system or a chamber loading error. | Reset the machine. Check:<br>a) that the door seal is clean and positioned correctly;<br>b) that there is water in the feed tank;<br>c) the amount of load in the chamber.<br>Then repeat the selected program.<br>If the problem persists, contact technical support. |

| Code   | Possible cause                                                                              | Possible solution                                                                                                                                                                                                                                                                 |
|--------|---------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ER 152 | There is a malfunction in the sterilisation condition holding/control system.               | Reset the machine. Check:<br>a) that there is water in the feed tank;<br>b) the amount of load in the chamber.<br>Then repeat the selected program.<br>If the problem persists, contact technical support.                                                                        |
| ER 153 | There is a malfunction in the sterilisation condition holding/control system.               |                                                                                                                                                                                                                                                                                   |
| ER 154 | There is a malfunction in the sterilisation condition holding/control system.               | Reset the machine. Check:<br>a) the tray support position and the tray positions in the chamber;<br>b) that there is water in the feed tank;<br>c) the amount of load in the chamber.<br>Then repeat the selected program.<br>If the problem persists, contact technical support. |
| ER 160 | There is a malfunction in the machine discharge system.                                     | Reset the machine. Check the amount of load in the chamber and then repeat the selected program.<br>If the problem persists, contact technical support.                                                                                                                           |
| ER 171 | There is a leak in the machine plumbing circuit or a malfunction in the vacuum pump.        | Reset the machine. Check:<br>a) that the door seal is clean and positioned correctly;<br>b) that the filter in the chamber is clean and positioned correctly.<br>Then repeat the selected program.<br>If the problem persists, contact technical support.                         |
| ER 172 | There is a malfunction in the machine drying system.                                        | Reset the machine.<br>Check the amount of load in the chamber and the Bio-X filter condition, then repeat the selected program.<br>If the problem persists, contact technical support.                                                                                            |
| ER 173 | There is a leak in the machine plumbing circuit or a malfunction in the vacuum pump.        | Reset the machine. Check:<br>a) that the door seal is clean and positioned correctly;<br>b) that the filter in the chamber is clean and positioned correctly.<br>Then repeat the selected program.<br>If the problem persists, contact technical support.                         |
| ER 180 | There is a malfunction in the pressure balancing system.                                    | Reset the machine.<br>Check the amount of load in the chamber and the Bio-X filter condition, then repeat the selected program.<br>If the problem persists, contact technical support.                                                                                            |
| ER 200 | There is a leak in the machine plumbing circuit or a malfunction in the vacuum pump.        | Reset the machine. Check:<br>a) that the door seal is clean and positioned correctly;<br>b) that the filter in the chamber is clean and positioned correctly.<br>Then repeat the selected program.<br>If the problem persists, contact technical support.                         |
| ER 201 | There is a leak in the machine plumbing circuit.                                            | Reset the machine.<br>Check that the door seal is clean and positioned correctly, then repeat the selected program.<br>If the problem persists, contact technical support.                                                                                                        |
| ER 202 | There is a leak in the machine plumbing circuit.                                            |                                                                                                                                                                                                                                                                                   |
| ER 203 | A Vacuum Test is running in chamber conditions other than ambient temperature.              | Reset the machine.<br>Repeat the Vacuum Test, making sure that the chamber is at ambient temperature. If the problem persists, contact technical support.                                                                                                                         |
| ER 204 | There is a malfunction in the pressure balancing system.                                    | Reset the machine.<br>Repeat the selected program.<br>If the problem persists, contact technical support.                                                                                                                                                                         |
| ER 310 | Leak in the hydraulic circuit of the machine or malfunctioning of the vacuum pump.          | Reset the machine. Check:<br>a) that the door seal is clean and positioned correctly;<br>b) that the filter in the chamber is clean and positioned correctly.<br>Then repeat the selected program.<br>If the problem persists, contact technical support.                         |
| ER 340 | Leak in the hydraulic circuit of the machine or malfunctioning of the water injection pump. | Reset the machine. Check:<br>a) that the door seal is clean and positioned correctly;<br>b) that there is water in the feed tank;<br>c) the amount of load in the chamber.<br>Then repeat the selected program.<br>If the problem persists, contact technical support.            |
| ER 341 | Leak in the hydraulic circuit of the machine or malfunctioning of the water injection pump. |                                                                                                                                                                                                                                                                                   |
| ER 342 | The safety valve has opened during test execution.                                          | Replace the safety valve.                                                                                                                                                                                                                                                         |
| ER 355 | The test execution time has been exceeded.                                                  | Check the condition of the safety valve and replace it if necessary.                                                                                                                                                                                                              |
| ER 360 | Malfunction in the machine discharge system.                                                | Reset the machine. Check the amount of load in the chamber and then repeat the selected program.<br>If the problem persists, contact technical support.                                                                                                                           |
| ER 370 | Malfunction in the machine drying system.                                                   | Reset the machine.<br>Check the amount of load in the chamber and the Bio-X filter condition. Then repeat the selected program.<br>If the problem persists, contact technical support.                                                                                            |
| ER 380 | Malfunction in the machine pressure balancing system.                                       |                                                                                                                                                                                                                                                                                   |
| AL 001 | The reset button was pressed.                                                               | Reset the machine.                                                                                                                                                                                                                                                                |
| AL 002 | A blackout has occurred.                                                                    | Reset the machine. Check:<br>a) that the power cable is connected correctly;<br>b) that there is mains power in the system.<br>Then repeat the selected program.<br>If the problem persists, contact technical support.                                                           |

| Code   | Possible cause                                                     | Possible solution                                                                                                                                                                                                                                                            |
|--------|--------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| AL 003 | There is a malfunction in the feed tank water filling system.      | Reset the machine. Check that the water filling hose is positioned and connected correctly, then repeat the selected procedure.<br>If the problem persists, contact technical support.<br>If water comes out of the overflow fitting at the rear, contact technical support. |
| AL 100 | There is a malfunction in the chamber probe reading circuit.       | Reset the machine.<br>Repeat the selected program.<br>If the problem persists, contact technical support.                                                                                                                                                                    |
| AL 101 | There is a malfunction in the generator probe reading circuit.     |                                                                                                                                                                                                                                                                              |
| AL 102 | There is a malfunction in the pipe bundle probe reading circuit.   |                                                                                                                                                                                                                                                                              |
| AL 103 | There is a malfunction in the pressure transducer reading circuit. | Reset the machine.<br>Repeat the selected program.<br>If the problem persists, contact technical support.                                                                                                                                                                    |
| AL 104 | There is a malfunction in the door closing/opening control system. |                                                                                                                                                                                                                                                                              |
| AL 105 | There is a malfunction in the door closing/opening control system. |                                                                                                                                                                                                                                                                              |
| AL 200 | There is a malfunction in the temperature control system.          |                                                                                                                                                                                                                                                                              |
| AL 201 | There is a malfunction in the temperature control system.          |                                                                                                                                                                                                                                                                              |
| AL 202 | There is a malfunction in the temperature control system.          | Reset the machine.<br>Contact technical support.                                                                                                                                                                                                                             |
| AL 203 | There is a malfunction in the pressure control system.             |                                                                                                                                                                                                                                                                              |
| AL 300 | There is a malfunction in the electronic system.                   | Reset the machine.<br>Repeat the selected program.<br>If the problem persists, contact technical support.                                                                                                                                                                    |
| AL 301 | There is a malfunction in the electronic system.                   |                                                                                                                                                                                                                                                                              |
| AL 302 | There is a malfunction in the electronic system.                   |                                                                                                                                                                                                                                                                              |
| AL 303 | There is a malfunction in the electronic system.                   |                                                                                                                                                                                                                                                                              |
| AL 304 | There is a malfunction in the electronic system.                   |                                                                                                                                                                                                                                                                              |
| AL 305 | There is a malfunction in the electronic system.                   |                                                                                                                                                                                                                                                                              |
| AL 306 | There is a malfunction in the electronic system.                   |                                                                                                                                                                                                                                                                              |
| AL 307 | There is a malfunction in the electronic system.                   |                                                                                                                                                                                                                                                                              |
| AL 308 | There is a malfunction in the electronic system.                   |                                                                                                                                                                                                                                                                              |
| AL 309 | There is a malfunction in the electronic system.                   |                                                                                                                                                                                                                                                                              |
| AL 310 | There is a malfunction in the electronic system.                   |                                                                                                                                                                                                                                                                              |
| AL 311 | There is a malfunction in the electronic system.                   |                                                                                                                                                                                                                                                                              |
| AL 401 | Generator malfunction.                                             |                                                                                                                                                                                                                                                                              |
| AL 402 | Band malfunction.                                                  |                                                                                                                                                                                                                                                                              |

## 6.6 Message solutions

| Code          | phase                            | Message                                        | Possible solution                                                                                           |
|---------------|----------------------------------|------------------------------------------------|-------------------------------------------------------------------------------------------------------------|
| <b>Msg001</b> | Standby or cycle start           | Feed tank empty                                | Fill the feed tank                                                                                          |
| <b>Msg002</b> | Standby or cycle start           | Feed tank full                                 | Empty the discharge tank                                                                                    |
| <b>Msg003</b> | Pin entry                        | Incorrect Pin                                  | Enter the correct Pin                                                                                       |
| <b>Msg004</b> | Water loading                    | Poor water quality                             | Change the water in the feed tank with better quality water                                                 |
| <b>Msg005</b> | Cycle start                      | Door open                                      | Close the door                                                                                              |
| <b>Msg008</b> | Vacuum Test start                | Chamber temperature over limit                 | Repeat Vacuum Test when the chamber temperature returns below the allowed limit (45°C)                      |
| <b>Msg009</b> | Machine switch-on                | Serial number inconsistency between CPU and MB | Contact The Dürr Dental Technical Service                                                                   |
| <b>Msg010</b> | Machine switch-on                | Lubricate the automated system of the door     | Lubricate the automated system according to the instructions in par. 5.1.9 of the user manual               |
| <b>Msg011</b> | Machine switch-on                | Replace Bio-X filter                           | Replace the Bio-X filter with original Dürr Dental spare part                                               |
| <b>Msg012</b> | Machine switch-on                | Conduct maintenance                            | Contact service technician for routine maintenance                                                          |
| <b>Msg013</b> | Machine switch-on                | Incorrect date                                 | Set clock/calendar of the machine again. If the problem persists, contact the Dürr Dental Technical Service |
| <b>Msg014</b> | Machine switch-on or cycle start | Do vacuum test                                 | Do vacuum test                                                                                              |
| <b>Msg016</b> | Pin entry                        | Pin too short                                  | Entry 4-digit PIN                                                                                           |
| <b>Msg017</b> | Pin entry                        | Duplicate PIN                                  | Entry correct PIN                                                                                           |
| <b>Msg018</b> | Machine switch-on                | MB FW revision incompatible with CPU FW        | Contact the Dürr Dental Technical Service                                                                   |

## 7 Technical specifications

|                                       | Hygoclave 50      |          |        | Hygoclave 50 Plus |          |        |
|---------------------------------------|-------------------|----------|--------|-------------------|----------|--------|
| Power supply voltage                  | 220 V             | 230 V    | 240 V  | 220 V             | 230 V    | 240 V  |
| Mains frequency                       | 60 HZ             | 50/60 HZ | 50 HZ  | 60 HZ             | 50/60 HZ | 50 HZ  |
| Power                                 | 1880 W            | 2050 W   | 2230 W | 2010 W            | 2200 W   | 2400 W |
| Maximum current consumption           | 8,5 A             | 8,9 A    | 9,3 A  | 9,1 A             | 9,6 A    | 10,0 A |
| Fuse (6.3x32)                         | T12 A H 250v      |          |        | T12 A H 250v      |          |        |
| Noise level                           | <60dB (A)         |          |        |                   |          |        |
| Electric shock protection class       | Class I           |          |        |                   |          |        |
| IP degree of protection               | IPX0              |          |        |                   |          |        |
| Operationg temperature                | From 5°C to 40°C  |          |        |                   |          |        |
| Chamber volume                        | approx. 17 litres |          |        | approx. 22 litres |          |        |
| Usable chamber volume                 | approx. 11 litres |          |        | approx. 14 litres |          |        |
| Feed tank volume                      | 5,2 litres        |          |        |                   |          |        |
| Weight per unit area of empty machine | 202 kN/m2         |          |        | 216 kN/m2         |          |        |
| Weight per unit area at full load     | 237 kN/m2         |          |        | 251 kN/m2         |          |        |
| Empty machine weight                  | 57 kg             |          |        | 61 kg             |          |        |
| Fully loaded machine weight           | 67 kg             |          |        | 71 kg             |          |        |
| Safety valve calibration pressure     | 2,5 bar           |          |        | 2,5 bar           |          |        |

### 7.1 Water specifications (distilled or demineralized)

|                      |                  |                    |                |
|----------------------|------------------|--------------------|----------------|
| Conductivity at 20°C | 0÷15 green µS/cm | 15÷40 yellow µS/cm | > 40 red µS/cm |
| Chlorides            | <2.0 mg/l        |                    |                |
| pH at 20°C           | 5.0 - 7.0        |                    |                |

## 8 Spare parts codes

6046100101 : Tray handle  
 6048100023 : 17 L tray  
 6048100024 : 22 L tray  
 6048100019 : Bio-X-Filter  
 6048100018 : Door seal  
 6048100015\*: Maintenance kit 1000 cycles  
 6048100138\*: Maintenance kit 5000 cycles 17 liters  
 6048100139\*: Maintenance kit 5000 cycles 22 liters

\* The annual maintenance must only be carried out by an authorised and qualified service technician

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## PT

- 

- 
- The screenshots show the 'Completo' screen in the '134° Standard' app. The top status bar shows '134° Standard', '11:20', and '15/03/2016'. The screen displays a thumbs-up icon, a clock showing '00:43:16', a pressure gauge showing '1.0 bar', and a temperature gauge showing '18 °C'. A green progress bar is at the bottom. The 'Completo' screen is divided into three sections: 'Lotto', 'Indici', and 'Carico'. Each section has a 'OK' button and a 'Fallito' button. The 'Fallito' button is highlighted with a red box in each section. Red circles with numbers 11, 14, and 16 indicate the sequence of actions.

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DUERR DENTAL FARO STERILIZATION S.r.l.