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China Council for the Promotion of International Trade  
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of GUANGDONG HARMONIC MEDICAL COMPANY LIMITED on the annexed DOCUMENT is genuine.

China Council for the Promotion  
of International Trade

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Dental Brushless Electric Motor

Operation Manual

Operations manuals for Dental Brushless Electric Motor:

MERCURY STAR 200EM

MERCURY STAR 400EM

MERCURY STAR 600EM

MERCURY STAR 800EM

Version No.: H02

Authorized representative on the territory of the Russian Federation –

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## Content

1. User instructions.....	5
2. Safety notes & Warnings .....	6
3. Appearance.....	7
4. Technical Specifications.....	7
5. Operation, transportation and storage conditions.....	8
6. Operation instructions.....	9
7. Maintenance.....	9
8. Marking.....	10
9. Materials in contact with human body.....	10
10. List of animal or human materials.....	10
11. Data on drugs and pharmaceutical substances in the medical unit.....	10
12. Complete set .....	10
13. Sterilization and disinfection .....	11
14. Service life .....	11
15. Warranty.....	11
16. Disposal and destruction procedure.....	12
17. Packing.....	12
18. List of international regulatory documents/ standards	Ошибка! Залкада не определена.
Annex 1. Installation scheme of model 200EM.....	12

Annex 2. EMC, radio frequencies requirements.....14  
Annex 2. The procedure for installing and checking the operation of the product..16  
**CONTACT DETAILS** .....23

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## 1. User instructions

### 1.1 User guide

Read the operation manual carefully prior to the first operation to avoid misuse and prevent from damage.

### 1.2 Target group

This unit is for dentists and dental qualified staff.

### 1.3 Symbol descriptions

The following symbols may appear in this manual, on the marking label, or on its accessories. Some of the symbols represent standards and compliances associated with the unit and its use.

CE Mark: conforms to essential requirements of the Medical Device Directive 93/42/EEC.

Date of manufacture.

Manufacturer

Only for indoor use

AC

Specifies serial number

Notes

Type B applied part

DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.

Refer to instruction manual / booklet

### 1.4 Classification

#### 1.4.1 Class II equipment;

#### 1.4.2 Type B applied part.

### 1.5. Purpose of the medical unit

The unit is intended to supply driving power to the dental handpiece for dental

treatment and operations in hospitals and dental clinics.

1.6. Operating conditions

It is used in medical institutions.

1.7. Field of use

Dentistry.

1.8. Potential user

This unit may only be used by dentists and qualified dental personnel familiar with the operation manual.

1.9. Indications

The unit is used to supply the driving power to the dental handpiece.

1.10. Contraindications

None.

1.11. Side effects

Not identified.

1.12. Risks of medical device application

During steps of risk analysis, it was concluded that this is a low-risk product and that any existing risks were eliminated through safety testing, proper selection of materials, and instructions for use. It was confirmed by several years of market experience with minimal complaints and no serious injuries. In addition, medical studies performed that during usage of the device did not detect serious injuries and showed a high level of success. Therefore, the product is safe.

## 2. Safety notes & Warnings

### 2.1 General information

The user must ensure that the unit works properly and is in satisfactory condition before each use or operation.

This unit is designed for use in dentistry only and may only be used by trained medical personnel. Any other type of use is not permitted.

### 2.2 Purpose-Proper use

“Proper use” includes following all the instructions for use and ensuring that all inspections and service tasks are performed.

2.2.1 Installation, expansions, adjustments, changes or repairs may only be carried out

by trained technicians. Before installation ensure that medical technician or distributor who offers the unit is qualified and authorized. In the case when all previously described actions were carried out by technician without a license and proper authorization, the manufacturer disclaims all the warranties.

2.2.2 When operating the unit, be sure to observe the instructions provided in the manual.

2.2.3 The unit must be cleaned and serviced according to instructions.

### 2.3 Warnings

2.3.1 Do not dispose of electrical appliances as unsorted municipal waste, use separate collection facilities. Contact your local government for information regarding the collection systems available. If electrical appliances were disposed in landfills or with domestic wastes, there is a risk that hazardous substances can leak into the groundwater and get into the food chain, damaging your health and well-being.

2.3.2 The device is not suitable for use in the presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide.

2.3.3 The user is prohibited to remove dental handpiece while motor is rotating.

2.3.4 It is prohibited to lubricate the micromotor with oil, avoid oil to enter inside of micromotor.

2.3.5 It is prohibited to disinfect micromotor by chemical disinfection.

## 3. Appearance


Fig.1 Control board

Fig.2 Micromotor and tubing

## 4. Technical Specifications

System: ----- Vector control system

Type: ----- Brushless motor

Revolutions: ----- 2,000 – 40,000 rpm  
 Torque: ----- 3.4Ncm  
 Power Input of adaptor: ----- AC 24 – 26.8 V  
 Dimension of motor: -----  $\phi 22 \times L71$ mm  
 Weight of motor: ----- 67 g  
 Noise: ----- <60dB  
 Lamp: ----- LED (white light)  
 Number of phases: ----- 3  
 Connection type: ----- ISO-E  
 Power: ----- 50 W  
 Altitude axis: ----- 52 mm  
 Axis type: ----- Asynchronous  
 Dimension of PCB: ----- L140 x W123 x H55mm  
 Total net weight: ----- 950.00 g

Micromotor tubing	
Length	Up to 1000 mm
Weight	200 g
Speed controller	
Dimensions	120×190 mm
Weight	150 g
Reverse Switch	
Dimensions	130×170 mm
Weight	250 g
Fastener	
Dimensions	120×110 mm
Weight	300 g

## 5. Operation, transportation and storage conditions

### 5.1 Conditions of transportation, storage and operation

Working environment temperature: -----  $-40^{\circ}\text{C} \sim +70^{\circ}\text{C}$

Relative humidity: ----- 10% ~ 100%

Atmosphere pressure: ----- 50Kpa ~ 106Kpa

## 6. Operation instructions

**Note:** Operation instructions for the electric motor is the same as operation instructions for air motor.

6.1 Pull out the motor from handpieces' holder of dental unit, press foot controller micro-motor starts and light indicator of motor lights on.

6.2 Press Speed controller to minus or plus the speed, speed down or up from 2000 to 40000 rpm.

6.3 Press Reverse switch to convert motor rotation direction.

**Warning:** Never attempt to convert micromotor rotation direction when micromotor is running.

## 7. Maintenance

7.1 Maintenance and repair may only be carried out by qualified and trained personnel.

7.2 Daily maintenance includes the micromotor cleaning with disinfecting alcohol-containing solution at the end of the working day, aqueous solution with up to 2,5% alcohol. Use seal caps at the front and back parts of micromotor before cleaning (as on the figure below).


**Warning:** Do not use perhydrol or similar disinfecting solution for cleaning. Do not oil or lubricate motor, it may cause motor failure and invalidate warranty replacement.

7.3 It is recommended to replace removable rubbers and ring gaskets after long-time use.

**Warning:** Do not sterilize motor by autoclaving.

**Warning:** It is forbidden to use chemical solutions for cleaning.

## 8. Marking

## 9. Materials in contact with human body

There is no direct contact with human body.

## 10. List of animal or human materials

There are no animal or human materials in the unit.

## 11. Data on drugs and pharmaceutical substances in the medical unit

There are no drugs and pharmaceutical substances in the unit.

## 12. Complete set

The delivery set includes:

Micromotor: ----- 1

Control board: ----- 1

Motor tubing: ----- 1

Speed controller: ----- 1

Reverse switch: ----- 1

Fasteners: ----- 1

User's manual: ----- 1

### 13. Sterilization and disinfection

The unit is non-sterile and does not require disinfection before use.

### 14. Service life

Service life for this product is 10 years.

### 15. Warranty

Manufacturer warrants its products to be in conformance with its own plans and specifications under normal use, transportation and storage.

**16. Disposal and destruction procedure**

The medical unit must be disposed in accordance with the requirements of the current legislation (including Sanitary-epidemiological requirements 2.1.7.2790-10).

**17. Packing**

Each component is individually packed in plastic bag and placed in cardboard box. Internal polyurethane inserts in transport package provides additional strength, prevents deformation and displacement while transportation.

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**18. List of regulatory documents / standards**

**International :**

EN ISO 13485: 2012	Medical Devices - Quality Management Systems - Regulatory Purpose Requirements
EN ISO 14971: 2012	Medical Devices I- Applying Risk Management to Medical Devices
EN 1041: 2008 + A1: 2013	Information provided by the manufacturer of medical devices
EN ISO 15223-1: 2012	Medical devices - Symbols to be used on the labels of medical devices, labeling, and information that must be provided a- Part 1: General requirements
EN 980: 2008	Designations for use in the labeling of medical devices
EN 60601-1: 2006 + A1: 2013	Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Basic Performance
EN 60601-1-2: 2007	"Medical electrical equipment - Part 1-2: General requirements for basic safety and essential operational quality m- Collateral Standard: Electromagnetic Compatibility i- requirements and tests"
EN 62304: 2006 + A1: 2015	Medical Device Software i- Software Life Cycle Processes
EN 60601-1-6: 2010 + A1: 2015	"Medical Electrical Equipment - Part 1-6: General Requirements for Basic Safety and Basic Operational Quality m- Supplementary Standard: Ease of Use"


**National**

No .	The code	Description
1	YY 0836 - 2011	Dental handpieces and micromotors
2	YY / T 1043 .1 -2016	Dentistry - integrated equipment
3	GB 9706.1 - 2007	Medical electrical equipment - Part 1: Basic safety requirements

4	YY 0505-2012	Medical electrical equipment - Part 1-2: Basic safety requirements and basic functional characteristics - Electromagnetic compatibility deadline - requirements and tests
5	GB / T 14710-2009	Basic requirements and test methods for medical electrical equipment
6	YY / T 0628-2008	Dentistry Device - Symbols Used for Labeling
7	GB / T4205-2010	Basic application and job security
8	GB5023.1-2008	Basic wiring requirements

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Annex 1. Installation scheme of model 200EM

**Attention:** for installation of the motor use the above scheme. The power should be connected as a last step.

According to the scheme the air valve need to be connected with T-valve for connection of water, air and spray.

**Note:** foot control is not a part of the product, foot control is a basic component for any dental unit. Please consult with supplier of dental units regarding the presence of the foot control.

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## Annex 2. EMC, radio frequencies requirements

<b>Manufacturer's recommendations and statement - electromagnetic radiation</b>		
The motor is intended for use in an electromagnetic environment as defined below. Customer or user of the motor must ensure its use in such environment		
<b>Exhaust emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment - recommendations</b>
RF emission CISPR 11	Group 1	The motor uses radio-frequency energy only for its internal functions. Therefore, its RF emissions are very low and cannot cause any interference to nearby electronic equipment  The motor is suitable for use in all institutions, including living premises and rooms directly connected to a public low voltage power network that supplies buildings used for household purposes.
RF emission CISPR 11	Class B	
emission, created by harmonic currents IEC 61000-3-2	No electric network input	
Voltage fluctuation/ power flickers IEC 61000-3-3	No electric network input	

<b>Manufacturer's recommendations and statement - susceptivness</b>			
The motor is intended for use in an electromagnetic environment as defined below. Customer or user of the motor must ensure its use in such environment			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance degree</b>	<b>Electromagnetic environment - recommendations</b>
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floor must be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic material, the relative humidity should be at least 30%.
Electric surges IEC 61000-4-4	± 2 kV for power supply lines  ± 1 kV for input/output lines	No electric network input  ± 1 kV for input/output lines	The quality of the electric power network should be typical for a commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line to line ± 2 kV line to line	No electric network input	The quality of the electric power network should be typical for a commercial or hospital environment.
Voltage dips, short interruptions and changes in supply voltage on the input power lines IEC 61000-4-11	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0.5 cycle  40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles  70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles	No electric network input	The quality of the electric power network should be typical for a commercial or hospital environment. If the OSSEOSTAP user requires continuous operation, it is recommended to power OSSEOSTAP from an uninterruptible power supply or battery.


	<5 % UT (>95 % dip in UT) for 5 sec		
Voltage frequency (50/60 Hz) of the magnetic field	3 A/m	30 A/m	The frequency of the magnetic field network should be at the level, characteristic of typical commercial or hospital environment
IEC 61000-4-8			
Note $U_T$ – AC mains voltage before applying the test level.			

Manufacturer's recommendations and statement - susceptivness			
The motor is intended for use in an electromagnetic environment as defined below. Customer or user of the motor must ensure its use in such environment			
Immunity test	IEC 60601 level tested	Compliance degree	Electromagnetic environment recommendations
Conducted RF/noise IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz	3 V	Portable and mobile radio frequency communications equipment should not be used near parts of the motor, including cables. The recommended distance is calculated by the equation applicable to the transmitter frequency.  <b>Recommended distance</b> $d = 1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz  $d = 2.3\sqrt{P}$ 800 MHz to 2.5GHz where p is the maximum rated output power of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended distance in meters (m). The field strength from fixed radio frequency transmitters, as determined by the investigation of the electromagnetic environment (a) at the site, should be less than the compliance level in each frequency range. (b) Interference may be present near equipment, marked with the following symbol: 
Note 1 At 80 MHz and 800 MHz, a higher frequency range is applied.			
Note 2 These recommendations may not apply in all situations. Absorption and reflection from structures, objects and people influence the electromagnetic waves propagation			
<sup>a</sup> The strength of the fixed transmitters field, such as base stations for radio (cellular / wireless) telephones and terrestrial mobile radio stations, amateur radio stations, AM and FM radio and television, can not be predicted theoretically accurately. To assess the electromagnetic environment based on fixed radio frequency transmitters, surveys of the electromagnetic situation on the site should be performed. If the measured field strength in the place, where the motor is used, exceeds the acceptable compliance level for RF Motor, which must be observed to verify normal operation. If faults are detected, additional measures			


such as reorienting or moving the motor are required  
b Field strength should be less than 3 V/m in the frequency range from 150 kHz to 80 MHz.

### Annex 3. The procedure for installing and checking the operation of the product.

Step 1: Open the control panel of the dental unit and locate it at appropriate place.

Step 2: Replace the standard hose for tool with the one for micromotor.

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

Step 3: Connect the water and air hoses of the micromotor to the corresponding hoses in the dental unit control panel using the air and water supply valve.


**Important:** The air and water supply valve is supplied separately of the micromotor.


Step 4: Connect the connector to the hose of micromotor on control board.


Step 5: Connect the control board to the power supply 24 V.

Step 6: Drill a small hole in the bottom of the dental unit control board to install the speed controller and reverse switch.


Step 7: Locate the micromotor in the tool holder of dental unit.

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
Step 8: Turn on the power of the dental unit and pull the micromotor out of the holder.

The micromotor is ready for use.

Step 9: Check the functions of the micromotor by pressing the pedal of the dental unit. When you press the pedal, the backlight of the micromotor should turn on and it should start rotating. Check the speed control and the reverse function using controllers that are fixed at the control panel. If any of the functions do not work, turn off the power of the unit and check correct connection of all wires. If everything is installed correctly and some functions or micromotor doesn't work, please contact the authorized representative of the manufacturer for warranty replacement of the

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## Content

1. User instructions.....	26
2. Safety notes & Warnings .....	27
3. Appearance.....	28
4. Technical Specifications .....	29
5. Operation, transportation and storage conditions.....	29
6. Operation instructions.....	30
7. Maintenance.....	30
8. Marking.....	31
9. Materials in contact with human body.....	31
10. List of animal or human materials.....	31
11. Data on drugs and pharmaceutical substances in the medical unit.....	31
12. Complete set .....	31
13. Sterilization and disinfection .....	32
15. Warranty.....	33
16. Disposal and destruction procedure.....	33
17. Packing.....	33
18. List of international regulatory documents/standards	<b>Ошибка! Закладка не определена.</b>
Annex 1. Installation scheme of model 400EM .....	37
Annex 2. Scheme of installation of the panel fastener .....	38
Annex 3. EMC, radio frequencies requirements.....	39

Annex 4. The procedure for installing and checking the operation of the product...41  
CONTACT DETAILS .....49

## 1. User instructions

### 1.1 User guide

Read the operation manual carefully prior to the first operation to avoid misuse and prevent from damage.

### 1.2 Target group

This unit is for dentists and dental qualified staff.

### 1.3 Symbol descriptions

The following symbols may appear in this manual, on the marking label, or on its accessories. Some of the symbols represent standards and compliances associated with the unit and its use.

CE Mark: conforms to essential requirements of the Medical Device Directive 93/42/EEC.

Date of manufacture.

Manufacturer

Only for indoor use

AC

Specifies serial number

Notes

Type B applied part

DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.

Refer to instruction manual / booklet

### 1.4 Classification

1.4.1 Class II equipment;

1.4.2 Type B applied part.

1.5. Purpose of the medical unit

The unit is intended to supply driving power to the dental handpiece for dental treatment and operations in hospitals and dental clinics.

1.6. Operating conditions

It is used in medical institutions.

1.7. Field of use

Dentistry.

1.8. Potential user

This unit may only be used by dentists and qualified dental personnel familiar with the operation manual.

1.9. Indications

The unit is used to supply the driving power to the dental handpiece.

1.10. Contraindications

None.

1.11. Side effects

Not identified.

1.12. Risks of medical device application

During steps of risk analysis, it was concluded that this is a low-risk product and that any existing risks were eliminated through safety testing, proper selection of materials, and instructions for use. It was confirmed by several years of market experience with minimal complaints and no serious injuries. In addition, medical studies performed that during usage of the device did not detect serious injuries and showed a high level of success. Therefore, the product is safe.

## 2. Safety notes & Warnings

### 2.1 General information

The user must ensure that the unit works properly and is in satisfactory condition before each use or operation.

This unit is designed for use in dentistry only and may only be used by trained medical personnel. Any other type of use is not permitted.

### 2.2 Purpose-Proper use

“Proper use” includes following all the instructions for use and ensuring that all

inspections and service tasks are performed.

2.2.1 Installation, expansions, adjustments, changes or repairs may only be carried out by trained technicians. Before installation ensure that medical technician or distributor who offers the unit is qualified and authorized. In the case when all previously described actions were carried out by technician without a license and proper authorization, the manufacturer disclaims all the warranties.

2.2.2 When operating the unit, be sure to observe the instructions provided in the manual.

2.2.3 The unit must be cleaned and serviced according to instructions.

### 2.3 Warnings

2.3.1 Do not dispose of electrical appliances as unsorted municipal waste, use separate collection facilities. Contact your local government for information regarding the collection systems available. If electrical appliances were disposed in landfills or with domestic wastes, there is a risk that hazardous substances can leak into the groundwater and get into the food chain, damaging your health and well-being.

2.3.2 The device is not suitable for use in the presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide.

2.3.3 The user is prohibited to remove dental handpiece while motor is rotating.

2.3.4 It is prohibited to lubricate the micromotor with oil, avoid oil to enter inside of micromotor.

2.3.5 It is prohibited to disinfect micromotor by chemical disinfection.

## 3. Appearance


Fig.1 Control board

Fig.2 Micromotor and tubing


Fig.3 Control panel

#### 4. Technical Specifications

System: -----	Vector control system
Type: -----	Brushless motor
Revolutions: -----	2,000 – 40,000 rpm
Torque: -----	3.4Ncm
Power Input of adaptor: -----	AC 24 – 26.8 V
Dimension of motor: -----	φ22×L71mm
Weight of motor: -----	67 g
Noise: -----	<60dB
Lamp: -----	LED (white light)
Number of phases: -----	3
Connection type:-----	ISO-E
Power: -----	60 W
Altitude axis: -----	46.3 mm
Axis type:-----	Asynchronous
Dimension of PCB:-----	L140 x W123 x H55mm
Total net weight: -----	1200.00 g

Micromotor tubing	
Length	Up to 1000 mm
Weight	200 g
Control panel	
Dimensions	300×200 mm
Weight	300 g
Panel fasteners	
Dimensions	150×180 mm
Weight	350 g
Fastener	
Dimensions	120×110 mm
Weight	300 g

## 5. Operation, transportation and storage conditions

### 5.1 Conditions of transportation, storage and operation

Working environment temperature: -----  $-40^{\circ}\text{C} \sim +70^{\circ}\text{C}$

Relative humidity: -----  $10\% \sim 100\%$

Atmosphere pressure: -----  $50\text{Kpa} \sim 106\text{Kpa}$

## 6. Operation instructions

**Note:** Operation instructions for the electric motor is the same as operation instructions for air motor.

6.1 Pull out the motor from handpieces' holder of dental unit, press foot controller micro-motor starts and light indicator of motor lights on.

6.2 Press buttons “+” or “-” on the control panel to adjust speed, it is possible to adjust it in range from 2000 to 40000 rpm, the speed is displayed on the panel.

6.3 Press Reverse button on the panel to convert motor rotation direction.

6.4 The panel has three buttons for automatic adjustment of the handpiece speed – 1:5, 1:1 and 1:20. When push the button information about the speed is displayed on the panel.

6.5 M1 and M2 buttons to save settings, when you push these buttons once again the motor will start to work according to settings.

**Warning:** Never attempt to convert micromotor rotation direction when micromotor is running.

## 7. Maintenance

7.1 Maintenance and repair may only be carried out by qualified and trained personnel.

7.2 Daily maintenance includes the micromotor cleaning with disinfecting alcohol-containing solution at the end of the working day, aqueous solution with up to 2,5% alcohol. Use seal caps at the front and back parts of micromotor before cleaning (as on the figure below).

**Warning:** Do not use perhydrol or similar disinfecting solution for cleaning. Do not oil or lubricate motor, it may cause motor failure and invalidate warranty replacement.

7.3 It is recommended to replace removable rubbers and ring gaskets after long-time use.

**Warning:** Do not sterilize motor by autoclaving.

**Warning:** It is forbidden to use chemical solutions for cleaning.

## 8. Marking

## 9. Materials in contact with human body

There is no direct contact with human body.

## 10. List of animal or human materials

There are no animal or human materials in the unit.

## 11. Data on drugs and pharmaceutical substances in the medical unit

There are no drugs and pharmaceutical substances in the unit.


## 12. Complete set

The delivery set includes:

Micromotor: -----|


Control board: -----|

Motor tubing: -----1

Control panel: -----1

Fasteners -----1

Panel fasteners: -----1

User's manual:-----1

### 13. Sterilization and disinfection

The unit is non-sterile and does not require disinfection before use.

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#### **14. Service life**

Service life for this product is 10 years.

#### **15. Warranty**

Manufacturer warrants its products to be in conformance with its own plans and specifications under normal use, transportation and storage.

#### **16. Disposal and destruction procedure**

The medical unit must be disposed in accordance with the requirements of the current legislation (including Sanitary-epidemiological requirements 2.1.7.2790-10).

#### **17. Packing**

Each component is individually packed in plastic bag and placed in cardboard box. Internal polyurethane inserts in transport package provides additional strength, prevents deformation and displacement while transportation.

## 18. List of regulatory documents / standards

### International :

EN ISO 13485: 2012	Medical Devices - Quality Management Systems - Regulatory Purpose Requirements
EN ISO 14971: 2012	Medical Devices I- Applying Risk Management to Medical Devices
EN 1041: 2008 + A1: 2013	Information provided by the manufacturer of medical devices
EN ISO 15223-1: 2012	Medical devices - Symbols to be used on the labels of medical devices, labeling, and information that must be provided a- Part I: General requirements
EN 980: 2008	Designations for use in the labeling of medical devices
EN 60601-1: 2006 + A1: 2013	Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Basic Performance
EN 60601-1-2: 2007	"Medical electrical equipment - Part 1-2: General requirements for basic safety and essential operational quality m- Collateral Standard: Electromagnetic Compatibility i- requirements and tests"
EN 62304: 2006 + A1: 2015	Medical Device Software i- Software Life Cycle Processes
EN 60601-1-6: 2010 + A1: 2015	"Medical Electrical Equipment - Part 1-6: General Requirements for Basic Safety and Basic Operational Quality m- Supplementary Standard: Ease of Use"

### National

No .	The code	Description
1	YY 0836 - 2011	Dental handpieces and micromotors
2	YY / T 1043 .1 -2016	Dentistry - integrated equipment
3	GB 9706.1 - 2007	Medical electrical equipment - Part 1: Basic safety requirements


4	YY 0505-2012	Medical electrical equipment - Part 1-2: Basic safety requirements and basic functional characteristics - Electromagnetic compatibility deadline - requirements and tests
5	GB / T 14710-2009	Basic requirements and test methods for medical electrical equipment
6	YY / T 0628-2008	Dentistry Device - Symbols Used for Labeling
7	GB / T 4205-2010	Basic application and job security
8	GB5023.1-2008	Basic wiring requirements

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[www.goszdravnadzor.ru](http://www.goszdravnadzor.ru)

Annex 1. Installation scheme of model 400EM


**Attention:** for installation of the motor use the above scheme. The power should be connected as a last step.

According to the scheme the air valve need to be connected with T-valve for connection of water, air and spray.

**Note:** foot control is not a part of the product, foot control is a basic component for any dental unit. Please consult with supplier of dental units regarding the presence of the foot control.

#### Annex 2. Scheme of installation of the panel fastener

### Annex 3. EMC, radio frequencies requirements

Manufacturer's recommendations and statement - electromagnetic radiation		
The motor is intended for use in an electromagnetic environment as defined below. Customer or user of the motor must ensure its use in such environment		
Exhaust emissions test	Compliance	Electromagnetic environment - recommendations
RF emission CISPR 11	Group 1	The motor is suitable for use in all institutions, including living premises and rooms directly connected to a public low voltage power network that supplies buildings used for household purposes.
RF emission CISPR 11	Class B	
emission, created by harmonic currents IEC 61000-3-2	No electric network input	
Voltage fluctuation/ power flickers IEC 61000-3-3	No electric network input	

Manufacturer's recommendations and statement - susceptivness			
The motor is intended for use in an electromagnetic environment as defined below. Customer or user of the motor must ensure its use in such environment			
Immunity test	IEC 60601 test level	Compliance degree	Electromagnetic environment - recommendations
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floor must be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic material, the relative humidity should be at least 30%.
Electric surges IEC 61000-4-4	± 2 kV for power supply lines  ± 1 kV for input/output lines	No electric network input  ± 1 kV for input/output lines	The quality of the electric power network should be typical for a commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line to line ± 2 kV line to line	No electric network input	The quality of the electric power network should be typical for a commercial or hospital environment.
Voltage dips, short interruptions and changes in supply voltage on the input power lines IEC 61000-4-11	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0.5 cycle  40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles  70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles	No electric network input	The quality of the electric power network should be typical for a commercial or hospital environment. If the OSSEOSTAP user requires continuous operation, it is recommended to power OSSEOSTAP from an uninterruptible power supply or battery.

	<5 % UT (>95 % dip in UT) for 5 sec		
Voltage frequency (50/60 Hz) of the magnetic field	3 A/m	30 A/m	The frequency of the magnetic field network should be at the level, characteristic of typical commercial or hospital environment
IEC 61000-4-8			
Note $U_T$ – AC mains voltage before applying the test level.			

**Manufacturer's recommendations and statement - susceptivness**

The motor is intended for use in an electromagnetic environment as defined below.  
Customer or user of the motor must ensure its use in such environment

Immunity test	IEC 60601 level tested	Compliance degree	Electromagnetic environment recommendations
Conducted RF/noise IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz	3 V	Portable and mobile radio frequency communications equipment should not be used near parts of the motor, including cables. The recommended distance is calculated by the equation applicable to the transmitter frequency.  <b>Recommended distance</b> $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5GHz where p is the maximum rated output power of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended distance in meters (m). The field strength from fixed radio frequency transmitters, as determined by the investigation of the electromagnetic environment (a) at the site, should be less than the compliance level in each frequency range. (b) Interference may be present near equipment, marked with the following symbol:
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	


Note 1 At 80 MHz and 800 MHz, a higher frequency range is applied.  
Note 2 These recommendations may not apply in all situations. Absorption and reflection from structures, objects and people influence the electromagnetic waves propagation

The strength of the fixed transmitters field, such as base stations for radio (cellular / wireless) telephones and terrestrial mobile radio stations, amateur radio stations, AM and FM radio and television, can not be predicted theoretically accurately. To assess the electromagnetic environment based on fixed radio frequency transmitters, surveys of the electromagnetic situation on the site should be performed. If the measured field strength in the place, where the motor is used, exceeds the acceptable compliance level for RF Motor, which must be observed to verify normal operation. If faults are detected, additional measures

such as reorienting or moving the motor are required  
b Field strength should be less than 3 V/m in the frequency range from 150 kHz to 80 MHz.


**Annex4 . The procedure for installing and checking the operation of the product.**

Step 1: Open the control panel of the dental unit and locate it at appropriate place.

Control board

Step 2: Replace the standard hose for tool with the one for micromotor.

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
[www.goszdravnadzor.ru](http://www.goszdravnadzor.ru)

Step 3: Connect the water and air hoses of the micromotor to the corresponding hoses in the dental unit control panel using the air and water supply valve.

**Important: The air and water supply valve is supplied separately of the micromotor.**


Valve


Step 4: Connect the connector to the hose of micromotor on control board.


Step 5: Connect the control board to the power supply 24 V.

Step 6: Drill a small hole to attach a remote control panel. Once you have been attached the special panel mount to install the control panel on it and connect it to the control board.

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

Step 7: Locate the micromotor in the tool holder of dental unit.


Step 8: Turn on the power of the dental unit and pull the micromotor out of the holder. The micromotor is ready for use.

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Step 9: Check the functions of the micromotor by pressing the pedal of the dental unit. When you press the pedal, the backlight of the micromotor should turn on and it should start rotating. Check the speed control and the reverse function using controllers that are fixed at the control panel. If any of the functions do not work, turn off the power of the unit and check correct connection of all wires. If everything is installed correctly and some functions or micromotor doesn't work, please contact the authorized representative of the manufacturer for warranty replacement of the equipment.

Pedal

## CONTACT DETAILS

### **Manufacturer of the device:**

#### **Guangdong Harmonic Medical Company Ltd**

Legal address: The Third Floor A, Building Nan, No. 9 of North Park East Road, Shishan Town, Nanhai District, Foshan City, Guangdong, China, 528225

Production address: The Third Floor A, Building Nan, No. 9 of North Park East Road, Shishan Town, Nanhai District, Foshan City, Guangdong, China, 528225

Phone: +86 757 85853030 Fax: +86 757 85853030

### **Authorized representative in the territory of the Russian Federation:**

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
## Content

1. User instructions.....	53
2. Safety notes & Warnings.....	54
3. Appearance.....	55
4. Technical Specifications.....	56
5. Installation.....	57
6. Operation instructions.....	58
7. Maintenance.....	59
8. Marking.....	59
9. Materials in contact with human body.....	60
10. List of animal or human materials.....	60
11. Data on drugs and pharmaceutical substances in the medical unit.....	60
12. Complete set.....	60
13. Sterilization and disinfection.....	61
14. Service life.....	61
15. Warranty.....	61
16. Disposal and destruction procedure.....	61
17. Packing.....	61
18. List of international regulatory documents/ standards	Ошибка! Залка не определена.
Annex 1. EMC, radio frequencies requirements.....	64

Annex 2. The procedure for installing and checking the operation of  
the product.....66  
CONTACT DETAILS ..... 72

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## 1. User instructions

### 1.1 User guide

Read the operation manual carefully prior to the first operation to avoid misuse and prevent from damage.

### 1.2 Target group

This unit is for dentists and dental qualified staff.

### 1.3 Symbol descriptions

The following symbols may appear in this manual, on the marking label, or on its accessories. Some of the symbols represent standards and compliances associated with the unit and its use.

	CE Mark: conforms to essential requirements of the Medical Device Directive 93/42/EEC.
	Date of manufacture.
	Manufacturer
	Only for indoor use
	AC
	Specifies serial number
	Notes
	Type B applied part
	DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.
	Refer to instruction manual / booklet

### 1.4 Classification

1.4.1 Class II equipment;

1.4.2 Type B applied part.

1.5. Purpose of the medical unit

The unit is intended to supply driving power to the dental handpiece for dental

treatment and operations in hospitals and dental clinics.

1.6. Operating conditions

It is used in medical institutions.

1.7. Field of use

Dentistry.

1.8. Potential user

This unit may only be used by dentists and qualified dental personnel familiar with the operation manual.

1.9. Indications

The unit is used to supply the driving power to the dental handpiece.

1.10. Contraindications

None.

1.11. Side effects

Not identified.

1.12. Risks of medical device application

During steps of risk analysis, it was concluded that this is a low-risk product and that any existing risks were eliminated through safety testing, proper selection of materials, and instructions for use. It was confirmed by several years of market experience with minimal complaints and no serious injuries. In addition, medical studies performed that during usage of the device did not detect serious injuries and showed a high level of success. Therefore, the product is safe.

## 2. Safety notes & Warnings

### 2.1 General information

The user must ensure that the unit works properly and is in satisfactory condition before each use or operation.

This unit is designed for use in dentistry only and may only be used by trained medical personnel. Any other type of use is not permitted.

### 2.2 Purpose-Proper use

“Proper use” includes following all the instructions for use and ensuring that all inspections and service tasks are performed.

#### 2.2.1 Installation, expansions, adjustments, changes or repairs may only be carried out

by trained technicians. Before installation ensure that medical technician or distributor who offers the unit is qualified and authorized. In the case when all previously described actions were carried out by technician without a license and proper authorization, the manufacturer disclaims all the warranties.

2.2.2 When operating the unit, be sure to observe the instructions provided in the manual.

2.2.3 The unit must be cleaned and serviced according to instructions.

### 2.3 Warnings

2.3.1 Do not dispose of electrical appliances as unsorted municipal waste, use separate collection facilities. Contact your local government for information regarding the collection systems available. If electrical appliances were disposed in landfills or with domestic wastes, there is a risk that hazardous substances can leak into the groundwater and get into the food chain, damaging your health and well-being.

2.3.2 The device is not suitable for use in the presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide.

2.3.3 The user is prohibited to remove dental handpiece while motor is rotating.

2.3.4 It is prohibited to lubricate the micromotor with oil, avoid oil to enter inside of micromotor.

2.3.5 It is prohibited to disinfect micromotor by chemical disinfection.

## 3. Appearance


Fig.1 Front view of Control unit

Fig.2 Micro-motor


Fig.3 Micromotor tubing


Fig.4 Rear view of control unit

**4. Technical Specifications**

- System: ----- Vector control system
- Type: ----- Brushless motor
- Revolutions: ----- 2,000 – 40,000 rpm
- Torque: ----- 3.4Ncm
- Power Input of adaptor: ----- AC 230V, 50Hz/60Hz
- Dimension of motor: -----  $\varnothing 22 \times L71$ mm
- Weight of motor: ----- 67 g
- Noise: ----- < 60dB
- Lamp: ----- LED (white light)
- Number of phases: ----- 3
- Power: ----- 60 W
- Altitude axis: ----- 46.3 mm
- Axis type: ----- Asynchronous
- Dimension of PCB: ----- L140 x W123 x H55mm
- Total net weight: ----- 950.00 g

Micromotor tubing	
Length	Up to 1000 mm
Weight	200 g
Adaptor	
Length	180×220 mm

Weight	200 g
Seal cap	
Dimensions	100×180 mm
Weight	300 g

## 5. Installation

### 5.1 Conditions of transportation, storage and operation

Working environment temperature:----- -40°C ~ +70°C

Relative humidity: ----- 10% ~ 100%

Atmosphere pressure: ----- 50Kpa ~ 106Kpa

### 5.2 Installation ways:

Step 1: Place Control unit on the top of Instrument tray of dental unit directly, and then connect 4 holes handpiece tube of dental unit to handpiece tube connector of Control unit.

Step 2: Connect Micromotor to micromotor tube with Control unit and placing Micromotor to handpiece holder of dental unit.

Step 3: Plug the adaptor to Control unit and connect to power.

Step 4: Power on the Control unit, taking out micromotor from holder, you can operate micromotor through Foot controller of dental unit. The operations ways are same as air motor.

**Note:** Foot controller is not included in complete set, Foot controller is a basic item for any dental unit. To ensure the availability of Foot controller consult the dental units' manufacturer.

## 6. Operation instructions

Marks and function of each button on the operation panel:


Fig.5 Operation panel

6.1 Press air pedal of Foot controller of Dental Unit, micro-motor starts and light indicator of motor lights on; Press water pedal of foot control of Dental Unit, the water supply to micromotor starts.

6.2 Press button Speed minus or Speed plus, speed down or up.

6.3 Press button REV/FWD to convert motor rotation direction. Indicator of R/F lights green when micromotor rotates forward; lights red when reverse.

**Warning:** Never attempt to convert micromotor rotation direction when micromotor is running.

6.4 There are 3 display models of rotation ratio, 1:5, 1:1 & 16:1. To switch rotation ratio, you need to press button Speed ratio to corresponding display model, corresponding Speed ratio indicator will light.

Rotation ratio	Speed range (rpm)
16:1	125~2,500
1:1	2,000~40,000
1:5	10,000~200,000

## 7. Maintenance

7.1 Maintenance and repair may only be carried out by qualified and trained personnel.

7.2 Daily maintenance includes the micromotor cleaning with disinfecting alcohol-containing solution at the end of the working day, aqueous solution with up to 2,5% alcohol. Use seal caps at the front and back parts of micromotor before cleaning (as on the figure below).

**Warning:** Do not use perhydrol or similar disinfecting solution for cleaning. Do not oil or lubricate motor, it may cause motor failure and invalidate warranty replacement.

7.3 It is recommended to replace removable rubbers and ring gaskets after long-time use.

**Warning:** Do not sterilize motor by autoclaving as it has negative impact on the motor lifetime and damage motor wiring.


Fig.6 Points of seal caps

## 8. Marking

**9. Materials in contact with human body**

There is no direct contact with human body.

**10. List of animal or human materials**

There are no animal or human materials in the unit.

**11. Data on drugs and pharmaceutical substances in the medical unit**

There are no drugs and pharmaceutical substances in the unit.

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
### Annex 1. EMC, radio frequencies requirements

Manufacturer's recommendations and statement - electromagnetic radiation		
The motor is intended for use in an electromagnetic environment as defined below. Customer or user of the motor must ensure its use in such environment		
Exhaust emissions test	Compliance	Electromagnetic environment - recommendations
RF emission CISPR 11	Group 1	The motor uses radio-frequency energy only for its internal functions. Therefore, its RF emissions are very low and cannot cause any interference to nearby electronic equipment  The motor is suitable for use in all institutions, including living premises and rooms directly connected to a public low voltage power network that supplies buildings used for household purposes.
RF emission CISPR 11	Class B	
emission, created by harmonic currents IEC 61000-3-2	No electric network input	
Voltage fluctuation/ power flickers IEC 61000-3-3	No electric network input	

Manufacturer's recommendations and statement - susceptibility			
The motor is intended for use in an electromagnetic environment as defined below. Customer or user of the motor must ensure its use in such environment			
Immunity test	IEC test level	Compliance degree	Electromagnetic environment - recommendations
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floor must be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic material, the relative humidity should be at least 30%.
Electric surges IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	No electric network input ± 1 kV for input/output lines	The quality of the electric power network should be typical for a commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line to line ± 2 kV line to line	No electric network input	The quality of the electric power network should be typical for a commercial or hospital environment.
Voltage dips, short interruptions and changes in supply voltage on the input power lines IEC 61000-4-11	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0.5 cycle  40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles  70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles	No electric network input	The quality of the electric power network should be typical for a commercial or hospital environment. If the OSSEOSTAP user requires continuous operation, it is recommended to power OSSEOSTAP from an uninterruptible power supply or battery.

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
<b>Manufacturer's recommendations and statement - susceptibility</b>			
The motor is intended for use in an electromagnetic environment as defined below. Customer or user of the motor must ensure its use in such environment			
<b>Immunity test</b>	<b>IEC 60601 level tested</b>	<b>Compliance degree</b>	<b>Electromagnetic environment recommendations</b>
Conducted RF/noise IEC 61000-4-6	3 $V_{rms}$ 150 kHz to 80 MHz	3 V	Portable and mobile radio frequency communications equipment should not be used near parts of the motor, including cables. The recommended distance is calculated by the equation applicable to the transmitter frequency.  <b>Recommended distance</b> $d = 1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz  $d = 2.3\sqrt{P}$ 800 MHz to 2.5GHz where p is the maximum rated output power of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended distance in meters (m). The field strength from fixed radio frequency transmitters, as determined by the investigation of the electromagnetic environment (a) at the site, should be less than the compliance level in each frequency range. (b) Interference may be present near equipment, marked with the following symbol:
Note 1 At 80 MHz and 800 MHz, a higher frequency range is applied.			
Note 2 These recommendations may not apply in all situations. Absorption and reflection from structures, objects and people influence the electromagnetic waves propagation			
<sup>a</sup> The strength of the fixed transmitters field, such as base stations for radio (cellular / wireless) telephones and terrestrial mobile radio stations, amateur radio stations, AM and FM radio and television, can not be predicted theoretically accurately. To assess the electromagnetic environment based on fixed radio frequency transmitters, surveys of the electromagnetic situation on the site should be performed. If the measured field strength in the place, where the motor is used, exceeds the acceptable compliance level for RF Motor, which must be observed to verify normal operation. If faults are detected, additional measures			

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
## Annex 2. The procedure for installing and checking the operation of the product.

Step 1: Unpack all components of the micromotor.


Step 2: Connect the power supply cable and micromotor cable to the micromotor unit.


Micromotor cable

Power supply cable

Step 3: Pull the tool hose that is already installed in the dental unit out of the micromotor hose socket where you want to locate the micromotor hose.

Step 4: Connect this tool hose to the connector on the micromotor unit.

Connector for tool hose


Step 5: Place the micromotor cable with the micromotor in the tool holder of which you pulled the tool hose out. Connect the power adapter to the network 220 V.

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

## CONTACT DETAILS

### **Manufacturer of the device:**

**Guangdong Harmonic Medical Company Ltd**

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
**MERCURY STAR 800EM**

Version No.: H02

Authorized representative on the territory of the Russian Federation:

LLC Novgorod

355003, Stavropol, Mira St., 367/21; +7(8652)525-888

[www.novgodent.ru](http://www.novgodent.ru)

## Contents

1. User instructions.....	75
2. Safety notes & Warnings .....	76
3. Appearance.....	78
4. Technical specifications .....	79
5. Installation .....	79
6. Operation Manual .....	81
7. Maintenance, repair, cleaning .....	83
8. Marking.....	84
9. Materials in contact with human body.....	84
10. List of animal or human materials.....	84
11. Data on the medicinal products and pharmaceutical substances contained in the medical device.....	84
12. Set.....	85
13. Sterilization, disinfection.....	86
14. Service life.....	86
15. Warranty.....	86
16. Disposal and destruction procedure.....	86
17. Packaging .....	86
18. List of international regulatory documents/ standards.....	86
Annex 1. EMC, radio frequencies requirements.....	90
CONTACT DETAILS.....	96

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	CE Mark: conforms to essential requirements of the Medical Device Directive 93/42/EEC.
	Date of manufacture
	Manufacturer
IPXO	General purpose equipment
IPXI	Splash protection type of the equipment
	Only for indoor use
	AC
SN	Specifies serial number
	Notes
	Type B applied part
	DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.
	Refer to the manual / booklet.

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2.2.1 Installation, expansions, adjustments, changes or repairs may only be carried out by trained technicians. Before installation ensure that medical technician or distributor who offers the unit is qualified and authorized. In the case when all

previously described actions were carried out by technician without a license and proper authorization, the manufacturer disclaims all the warranties.

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
Fig.4 Control unit, rear view

#### 4. Technical specifications

System: ----- Vector control system  
Type: ----- Brushless motor  
Rotation speed: ----- 2,000 – 40,000 rpm  
Torque: ----- 3.4 Ncm  
Input power: ----- AC 230V, 50/60 Hz  
Motor dimensions: -----  $\varnothing 22 \times L71$  mm  
Motor weight: ----- 67 gr.  
Noise: ----- < 60 dB  
Light source: ----- LED (white light)  
Number of phases: ----- 3  
Power: ----- 60 W  
Rotation axis height: ----- 46.3 mm  
Axis type: ----- Asynchronous  
Operation unit size: ----- L150.00 x W140.00 x H60.00 mm  
Total net weight: ----- 1250.00 gr.

Micromotor tubing	
Length	Up to 1000 mm
Weight	200 g
Adaptor	
Length	180×220 mm
Weight	200 g
Seal cap	
Dimensions	100×180 mm
Weight	300

#### 5. Installation

##### 5.1 Conditions of transportation, storage and operation

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
## 6. Operation Manual

Marking ad functions of each button on the working panel:

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**Optional:** You can use a micromotor without a dental unit, the micromotor can work independently. Press the power button on the working panel, the micromotor will start and the indicator will light up, then proceed according to items 6.2, 6.3, 6.4.

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## 10. List of animal or human materials

The device does not contain materials of animal and (or) human origin.

## 11. Data on the medicinal products and pharmaceutical substances contained in the medical device

The device does not contain medicinal products and pharmaceutical substances.

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

**14. Service life**

The service life of the product is 10 years.

**15. Warranty**

The manufacturer guarantees the conformity of the product to the declared characteristics, provided that the requirements of operation, transportation and storage are observed.

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The medical device must be disposed of in accordance with the requirements of the current legislation (including SanPiN 2.1.7.2790-10).

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**18. List of international regulatory documents/ standards**

EN ISO 13485: 2012	Medical Devices - Quality Management Systems - Requirements For Regulatory Purposes
EN ISO 14971: 2012	Medical Products - Application Of Risk Management To Medical Devices
EN 1041:2008+A1: 2013	Information Provided By The Manufacturer Of Medical Devices
EN ISO 15223-1:2012	Medical Products - Symbols To Be Used With Medical Device Labels, Labelling And Information To Be Supplied -- Part 1: General requirements
EN 980:2008	Symbols For Use In The Labeling Of Medical Devices
EN 60601-1:2006+A1:2013	Medical Electrical Equipment - Part 1: General Requirements For Basic Safety And Essential Performance.
EN 60601-1-2:2007	Medical Electrical Equipment - Part 1-2: General Requirements For Basic Safety And Essential Performance - Collateral standard: Electromagnetic compatibility - Requirements and tests
EN 62304:2006+A1:2015	Medical device software. Software life-cycle processes
EN 60601-1-6:2010+A1:2015	Medical Electrical Equipment - Part 1-6: General requirements for basic safety and essential performance. Collateral standard: Usability

### Annex 1. EMC, radio frequencies requirements.

Manufacturer's recommendations and statement - electromagnetic radiation			
The motor is intended for use in an electromagnetic environment as defined below. Customer or user of the motor must ensure its use in such environment			
Exhaust test	emissions	Compliance	Electromagnetic environment - recommendations
RF CISPR 11	emission	Group 1	The motor uses radio-frequency energy only for its internal functions. Therefore, its RF emissions are very low and cannot cause any interference to nearby electronic equipment  The motor is suitable for use in all institutions, including living premises and rooms directly connected to a public low voltage power network that supplies buildings used for household purposes.
RF CISPR 11	emission	Class B	
emission, created by harmonic currents		No electric network input	
Voltage fluctuation/power flickers		No electric network input	
IEC 61000-3-2			
IEC 61000-3-3			


Manufacturer's recommendations and statement - susceptivness			
The motor is intended for use in an electromagnetic environment as defined below. Customer or user of the motor must ensure its use in such environment			
Immunity test	IEC test level	Compliance degree	Electromagnetic environment - recommendations
Electrostatic discharge (ESD)	$\pm 6$ kV contact $\pm 8$ kV air	$\pm 6$ kV contact $\pm 8$ kV air	Floor must be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic material, the relative humidity should be at least 30%.
IEC 61000-4-2			
Electric surges	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/output lines	No electric network input $\pm 1$ kV for input/output lines	The quality of the electric power network should be typical for a commercial or hospital environment.
IEC 61004-4			
Surge	$\pm 1$ kV line to line $\pm 2$ kV line to line	No electric network input	The quality of the electric power network should be typical for a commercial or hospital environment.
IEC 61000-4-5			
Voltage dips, short interruptions and changes in supply voltage on the input lines	$<5$ % $U_T$ ( $>95$ % dip in $U_T$ ) for 0.5 cycle  $40$ % $U_T$ ( $60$ % dip in $U_T$ ) for 5 cycles  $70$ % $U_T$ ( $30$ % dip in $U_T$ ) for 25 cycles  $<5$ % $U_T$ ( $>95$ % dip in $U_T$ ) for 5 sec	No electric network input	The quality of the electric power network should be typical for a commercial or hospital environment. If the OSSEOSTAP user requires continuous operation, it is recommended to power OSSEOSTAP from an uninterruptible power supply or battery.
IEC 61000-4-11			

Voltage frequency (50/60 Hz) of the magnetic field	3 A/m	30 A/m	The frequency of the magnetic field network should be at the level, characteristic of typical commercial or hospital environment
IEC 61000-4-8			
Note $U_T$ – AC mains voltage before applying the test level.			


Manufacturer's recommendations and statement - susceptibility			
The motor is intended for use in an electromagnetic environment as defined below. Customer or user of the motor must ensure its use in such environment			
Immunity test	IEC 60601 level tested	Compliance degree	Electromagnetic environment recommendations
Conducted RF/noise IEC 61000-4-6	3 $V_{rms}$ 150 kHz to 80 MHz	3 V	Portable and mobile radio frequency communications equipment should not be used near parts of the motor, including cables. The recommended distance is calculated by the equation applicable to the transmitter frequency.  <b>Recommended distance</b> $d = 1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz  $d = 2.3\sqrt{P}$ 800 MHz to 2.5GHz where p is the maximum rated output power of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended distance in meters (m). The field strength from fixed radio frequency transmitters, as determined by the investigation of the electromagnetic environment (a) at the site, should be less than the compliance level in each frequency range. (b) Interference may be present near equipment, marked with the following symbol:
Note 1 At 80 MHz and 800 MHz, a higher frequency range is applied.			
Note 2 These recommendations may not apply in all situations. Absorption and reflection from structures, objects and people influence the electromagnetic waves propagation			
<sup>a</sup> The strength of the fixed transmitters field, such as base stations for radio (cellular / wireless) telephones and terrestrial mobile radio stations, amateur radio stations, AM and FM radio and television, can not be predicted theoretically accurately. To assess the electromagnetic environment based on fixed radio frequency transmitters, surveys of the electromagnetic situation on the site should be performed. If the measured field strength in the place, where the motor is used, exceeds the acceptable compliance level for RF Motor, which must be observed to verify normal operation. If faults are detected, additional measures such as reorienting or moving the motor are required <sup>b</sup> Field strength should be less than 3 V/m in the frequency range from 150 kHz to 80 MHz.			

**Appendix 2. The procedure for installing and checking the operation of the product.**


Step 1: Unpack all components of the micromotor.


Step 2: Connect the power supply cable and micromotor cable to the micromotor unit.

Step 3: Pull the tool hose that is already installed in the dental unit out of the micromotor hose socket where you want to locate the micromotor hose.

Step 4: Connect this tool hose to the connector on the micromotor unit.

Connector for tool hose

Step 5: Place the micromotor cable with the micromotor in the tool holder of which you pulled the tool hose out. Connect the power adapter to the network 220 V.

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


#### **CONTACT DETAILS**

##### **Manufacturer of the device:**

**Guangdong Harmonic Medical Company Ltd**

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Phone: +86 757 85853030      Fax: +86 757 85853030

##### **Authorized representative in the territory of the Russian Federation:**

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